



LC TECHNOLOGY SOLUTIONS INC.

2C Fanaras Drive, Salisbury, MA 01952 • Phone: 978-255-1620 • Fax: 978-428-0222

E-mail: info@lctechinc.com • Web: www.lctechinc.com

LC-1 Stand Alone Glovebox with Touch Screen

Operation/Installation Manual





Table of Contents

Section 1: System Overview	5
Section 2: Caution / Warning Information	7
Section 3: System Arrival	8
3.1 Disassembly of Crate	8
3.2 Remove Loose Parts	9
3.3 Remove System from Crate	10
3.4 Remove Gas Purifier Cart	11
Section 4: Assembly	12
4.1 Matching Labels	12
Section 5: Hook Up Gas / Power	13
5.1 Gas Connections	13
5.1 Power Connections	14
Section 6: Venting	15
Section 7: Window Removal	16
Section 8: Window Replacement	17
Section 10: Changing a Glove without a Glove Port Cover	19
Section 11: Leak Checking	23
Section 12: Purging	24
Section 13: Operational Instructions	25
13.1 Main Screen	25
13.2 Pressure Control	26
13.3 Automatic Purge	27
13.4 Automatic Purge Function Connected to Oxygen Level	30
13.5 Circulation	32
Section 14: Antechamber Operation	33
14.1 Bringing Items into the Glovebox	33
14.2 Removing Items from Glovebox	38
14.3 Automatic Antechamber Control / Loadlock	39
14.4 Manual Purge for Large Antechamber	45
14.5 Mini Antechamber Operation	46



LC TECHNOLOGY SOLUTIONS INC.

2C Fanaras Drive, Salisbury, MA 01952 • Phone: 978-255-1620 • Fax: 978-428-0222
E-mail: info@lctechinc.com • Web: www.lctechinc.com

14.6 Mini Antechamber with Automatic Control for Door Lock.....	48
14.7 Manual Purge for Mini Antechamber	50
Section 15: Regeneration Mode	51
15.1 Testing the Regeneration Process	55
15.2 Regeneration of a Dual Filter Column System	56
Section 16: Chamber Settings.....	57
16.1 Set Points	57
16.2 Enable / Disable.....	60
16.3 Dual Filter Columns	61
Section 17: Analyzers	62
Section 18: Solvent Removal Systems, Operation and Maintenance	63
18.1 Manual Solvent Removal System Operation	63
18.2 Manual Solvent Removal System Maintenance	64
18.3 Automatic Solvent Removal System Operation.....	66
18.4 Automatic Solvent Reactivation.....	67
Section 19: Internal Charcoal Trap.....	70
Section 20: Freezer Operation and Maintenance	72
20.1 Freezer Operation.....	72
20.2 Freezer Maintenance	73
Section 21: Alarm Messages.....	74
Section 22: Manual Purge	78
Section 23: Maintenance Schedule & Recommended Spare Parts	80
Section 24: Operational Accessories.....	81
24.1 Laminar Flow	81
24.2 UV Ozone Cleaner	83
24.3 Precision Hot Plate	86
24.4 Spin Coater.....	88
24.5 UV Press.....	89
24.6 FUJIFILM Dimatix Materials Printer.....	94
24.7 Fisnar Glue Dispensing Robot	95
24.8 Heat & Cooling Plate.....	96
24.9 Heated Tray	98
24.9a Heated Tray Controller Set-Up.....	100
24.10 Solvent Sensor.....	102



LC TECHNOLOGY SOLUTIONS INC.

2C Fanaras Drive, Salisbury, MA 01952 • Phone: 978-255-1620 • Fax: 978-428-0222
E-mail: info@lctechinc.com • Web: www.lctechinc.com

24.11 Power Saver Mode	104
Section 25: Technical Support Contact Information from LC Technology Solutions & Partners	105




LC TECHNOLOGY SOLUTIONS INC.

2C Fanaras Drive, Salisbury, MA 01952 • Phone: 978-255-1620 • Fax: 978-428-0222
E-mail: info@lctechinc.com • Web: www.lctechinc.com

LC-1 Stand Alone Glovebox with Touch Screen Operation Manual

Section 1: System Overview




Gas Purification System



LC TECHNOLOGY SOLUTIONS INC.

2C Fanaras Drive, Salisbury, MA 01952 • Phone: 978-255-1620 • Fax: 978-428-0222

E-mail: info@lctechinc.com • Web: www.lctechinc.com



The Project Number is located on a sticker on the back of your system. Please have this number available when calling for service information.



Section 2: Caution / Warning Information

1. Glovebox system is heavy. Use a fork lift to remove from pallet upon arrival, if available.
2. Make sure inert gas is regulated at 60-80 psi.
3. Make sure regeneration gas is regulated at 15 psi.
4. Common vent line should be vented.
5. It is very important to not run out of gas during purging of the glovebox.
6. Do not over tighten antechamber door handles.
7. Regeneration cycle times are critical to the system running properly. Do not change these settings.
8. Over / under pressure alarms in the system are dangerous and caution should be taken not to reach them.
9. Do not open glovebox with inert gas on.
10. Low oxygen hazard inside glovebox.
11. Make sure glovebox is well ventilated before entering.



LC TECHNOLOGY SOLUTIONS INC.

2C Fanaras Drive, Salisbury, MA 01952 • Phone: 978-255-1620 • Fax: 978-428-0222
E-mail: info@lctechinc.com • Web: www.lctechinc.com

Section 3: System Arrival

3.1 Disassembly of Crate



1. System will arrive in a crate.
2. Remove crate panels in the following order.
 - a. Top panel
 - b. Shipping braces, if applicable
 - c. Large left & right panels
 - d. Small front & back panels




3. Use a cordless drill with a Phillips head driver bit to remove screws.



4. This is how the system will look inside the crate.



3.2 Remove Loose Parts



1. Cut banding straps and remove boxes containing loose parts for glovebox system.




2. Parts will come packed in bubble wrap. Remove bubble wrap.



LC TECHNOLOGY SOLUTIONS INC.

2C Fanaras Drive, Salisbury, MA 01952 • Phone: 978-255-1620 • Fax: 978-428-0222
E-mail: info@lctechinc.com • Web: www.lctechinc.com


3.3 Remove System from Crate



1. Once sides are off, unbolt lag bolts from Z-brackets using a $\frac{3}{4}$ " wrench.

OR

Remove shipping braces with phillips head screw driver.




2. Lower stand onto its wheels by adjusting leveling feet with wrench.



3. Remove system from crate.

CAUTION: System is heavy. Use a fork lift to remove from pallet, if available.




4. Carefully remove shrink wrap.



LC TECHNOLOGY SOLUTIONS INC.

2C Fanaras Drive, Salisbury, MA 01952 • Phone: 978-255-1620 • Fax: 978-428-0222
E-mail: info@lctechinc.com • Web: www.lctechinc.com

3.4 Remove Gas Purifier Cart



1. Remove/cut banding straps.
2. Unscrew 2x4's holding purifier in place. Use a cordless drill with a phillips head driver bit to remove screws.



3. Remove purifier from crate. Use a fork lift to remove from pallet, if available.



4. Carefully remove shrink wrap.



Section 4: Assembly

4.1 Matching Labels



1. Various components of the glovebox will need to be attached to the system. These components have been number labeled to show connection points. Match corresponding numbers (i.e. 1 to 1, 2 to 2, and so on).


See examples to the left.



Section 5: Hook Up Gas / Power


5.1 Gas Connections

System will be supplied with reinforced 3/8" Tygon tubing already attached to the gas purifier.



1. Find the tube labeled Gas and connect this tube to your inert gas (nitrogen or argon) supply.

NOTE: The 3/8" Tygon tube should be connected to inert gas at 60-80 psi.



2. Find the tube labeled Regas and connect this tube to your regeneration gas supply.

NOTE: The 3/8" Tygon tube should be connected to regeneration gas at 15 psi.

Regeneration gas consists of 3-5% hydrogen, balance is nitrogen or argon.

Regeneration gas flow rate is 15 liters per minute.



LC TECHNOLOGY SOLUTIONS INC.

2C Fanaras Drive, Salisbury, MA 01952 • Phone: 978-255-1620 • Fax: 978-428-0222
E-mail: info@lctechinc.com • Web: www.lctechinc.com

5.1 Power Connections

1. System will have 115V/220V power cord coming from the back of the purifier.
2. Plug cord into a standard 115V, 15 Amp outlet or standard 220V outlet for international use.



Electrical Feedthroughs:

US locations plug electrical feedthrough into 115V power supply.

International locations plug electrical feedthrough into 220V power supply.



Vacuum Pump:

LC-1 Gloveboxes with an RGP-1 Gas Purifier come prewired into the glovebox system.



LC TECHNOLOGY SOLUTIONS INC.

2C Fanaras Drive, Salisbury, MA 01952 • Phone: 978-255-1620 • Fax: 978-428-0222
E-mail: info@lctechinc.com • Web: www.lctechinc.com

Section 6: Venting



Common Vent Line:

Systems will come with a common vent line. All items on the system will be plumbed to this common vent line. This vent line will have a 1 ½" OD diameter connection point and exhaust 8 CFM.


NOTE: It is highly recommended that this line be vented to an exhaust system.

Vent Line Port



Preferred Method of Vent Connection:

Vent the glovebox using a loose fit duct work adapter (Snorkel Connection).



NOTE: If System was supplied with an over pressure relief device. Do not place the fume hood connection directly on the automatic pressure relief / exhaust valve. Place the connection over the valve as shown.

Placing the connection directly on top of the automatic pressure relief valve will cause it to malfunction.



LC TECHNOLOGY SOLUTIONS INC.


2C Fanaras Drive, Salisbury, MA 01952 • Phone: 978-255-1620 • Fax: 978-428-0222
E-mail: info@lctechinc.com • Web: www.lctechinc.com

Section 7: Window Removal


NOTE: The window can be removed for placing equipment in the glovebox that is too big to be brought in through the antechamber.



1. Undo all star knobs with the exception of two (2) at the top and two (2) at the bottom.



2. Remove the remaining star knobs being careful not to let the window and window frame fall.



3. Remove the window frame and set aside.




4. Remove window and complete required work.




Section 8: Window Replacement




1. Place window on bottom window studs and push window forward into position against gasket.



2. Replace window frame.



3. Replace top two (2) and bottom two (2) star knobs.



4. Replace remaining star knobs.

NOTE: Do not tighten until all knobs have been started.



5. Tighten star knobs until window frame contacts glovebox.




LC TECHNOLOGY SOLUTIONS INC.

2C Fanaras Drive, Salisbury, MA 01952 • Phone: 978-255-1620 • Fax: 978-428-0222
E-mail: info@lctechinc.com • Web: www.lctechinc.com


Section 9: Glove Attachment




1. Place glove onto glove port.



2. Place hand in glove and align glove hand in a comfortable position with thumbs facing in an upward position.
3. Move the glove forward until it meets the inner most glove port groove.



4. Place first glove O-ring in center glove port groove.



5. Place second glove O-ring in outer most glove port groove.



LC TECHNOLOGY SOLUTIONS INC.

2C Fanaras Drive, Salisbury, MA 01952 • Phone: 978-255-1620 • Fax: 978-428-0222
E-mail: info@lctechinc.com • Web: www.lctechinc.com

Section 10: Changing a Glove without a Glove Port Cover






LC TECHNOLOGY SOLUTIONS INC.


2C Fanaras Drive, Salisbury, MA 01952 • Phone: 978-255-1620 • Fax: 978-428-0222
E-mail: info@lctechinc.com • Web: www.lctechinc.com



1. Remove O-ring nearest to window.



2. Glove with O-ring removed.




3. Fold glove back over existing O-ring as shown to the left.




LC TECHNOLOGY SOLUTIONS INC.


2C Fanaras Drive, Salisbury, MA 01952 • Phone: 978-255-1620 • Fax: 978-428-0222
E-mail: info@lctechinc.com • Web: www.lctechinc.com



4. Place new glove over the old glove. Make sure the new glove is mounted in O-ring groove.



5. Replace O-ring nearest to window.




6. Reach into new glove and carefully grab folded back edge of old glove.




LC TECHNOLOGY SOLUTIONS INC.

2C Fanaras Drive, Salisbury, MA 01952 • Phone: 978-255-1620 • Fax: 978-428-0222
E-mail: info@lctechinc.com • Web: www.lctechinc.com



7. Pull folded back edge of old glove until the old glove and O-ring fall into the glovebox.



8. Reattach second O-ring and glove change is complete.

NOTE: It is recommended to have circulation off during this procedure. You will need to purge the system for ten (10) minutes after the glove change is complete. Once purge is complete circulation can be turned back on.



Section 11: Leak Checking

After assembly of glovebox system it is important to leak check the system prior to purging.

NOTE: Leak checking is an important step in making sure your glovebox functions properly.

1. Turn on power to system.
2. To change settings refer to Section 16.1 Set Points.
3. Change Chamber Pressure Control settings to the following set points: 6.5 mBAR high and 2.0 mBAR low set points. These set points are for testing purposes only.
4. Press Return to return to Main Menu. Press Chamber Pressure Control button. ON should be displayed.
5. Gas will begin to flow into the box until the low set point is reached.
6. Using the foot pedal, press gas pedal to increase the pressure to the upper set point.

NOTE: When the upper set point is reached the vacuum pump will turn on lowering the pressure. Let the glovebox pressure settle for a couple of minutes.

If you are leak checking your glovebox after your system has been running, turn OFF Circulation for 30 Minutes prior to going on to Step #7 below.

7. After the glovebox pressure has settled, time how fast pressure drops 0.1 mBAR. For a successful test the pressure should not drop more than 0.1 mBAR for three (3) minutes.

e.g.: When the glovebox is at 2.5 mBAR it should hold between 2.4 and 2.5 mBAR for three (3) minutes.

NOTE: If the test does not reach the above requirement go through the fittings to check for a leak or any loose fittings. Once fittings have been checked repeat test until successful.

8. After successful leak check return pressure settings to 2.5 mBAR high and 0.5 mBAR low set points.



Section 12: Purging

After successfully leak checking the glovebox it is very important to purge the room air from the system. This will ensure your system works properly.

The chart below shows the recommended purge times and gas needed for each glovebox model at start-up.

Time and Gas Required:		
Glovebox Model No.: LCPW-125	0.5 cylinder, 150 c/f	20 minutes
Glovebox Model No.: LCPW-225	1.0 cylinder, 300 c/f	40 minutes
Glovebox Model No.: LC-100	2.5 cylinder, 750 c/f	2 hours
Glovebox Model No.: LC-150	3.5 cylinder, 1050 c/f	2 hours 30 min
Glovebox Model No.: LC-180	4.0 cylinder, 1200 c/f	3 hours
Glovebox Model No.: LC-200	5.0 cylinder, 1500 c/f	3 hours 30 min

Refer to Section 13.3 Automatic Purge and Section 21 Manual Purge for detailed instructions on how to purge the glovebox.

NOTE: 240 Minutes is the maximum amount of time you can set auto purge to run.

NOTE: You should only set auto purge to 240 minutes if you are using a house supply or dewar of inert gas when purging.

NOTE: You should only set auto purge to 45 minutes if you are using cylinders of inert gas when purging.

NOTE: Nitrogen/Argon has to be 99.995% or better.

CAUTION: It is very important to not run out of gas during purging of the glovebox.

NOTE: Once you have successfully completed purging your glovebox you can then turn circulation ON. (Please see Section 13.4, page 29) The oxygen and moisture levels displayed are only accurate when circulation is ON. The analyzers are flow sensitive and the circulation blower is used to flow gas through the analyzers to get a sample of gas from the glovebox. When first turning circulation on depending on the size of your glovebox, oxygen and moisture levels can take several hours to a day to read low levels. During this time they should be trending downward.




Section 13: Operational Instructions


13.1 Main Screen

Chamber Pressure Control ON is the everyday operating mode for the glovebox. Chamber Pressure Control starts the pressure control function and allows for glovebox circulation.

When the system is turned on, the display will show the following screen:



This is the main screen. All glovebox functions can be accessed through this screen.



Press Chamber Pressure Control to enable pressure control and glovebox circulation.

ON will be displayed above the Chamber Pressure Control button.



13.2 Pressure Control

1. From the Main Screen enable Chamber Pressure Control.
2. To enable press Chamber Pressure Control button. ON will be displayed above the Chamber Pressure Control button.
3. Pressure settings have been preset at the factory to a low limit of 0.5 mBAR and a high limit of 2.5 mBAR.

NOTE: This is a typical pressure range for the glovebox.

4. The system has been supplied with foot pedals to help control the pressure between the high and low set points.

Use the foot pedals to help control the pressure when inserting and removing hands from gloves.




LC TECHNOLOGY SOLUTIONS INC.

2C Fanaras Drive, Salisbury, MA 01952 • Phone: 978-255-1620 • Fax: 978-428-0222

E-mail: info@lctechinc.com • Web: www.lctechinc.com

13.3 Automatic Purge

NOTE: For Manual Purge instructions refer to Section 22.



LC TECHNOLOGY SOLUTIONS EL103

CHAMBER PRESSURE: +0000.0 mBAR

CHAMBER OXYGEN LEVEL: +0000.0 PPM

CHAMBER MOISTURE LEVEL: +0000.0 PPM

OFF

CHAMBER LIGHTS

OFF

CHAMBER PRESSURE CONTROL

CHAMBER CIRCULATION

CHAMBER SETTINGS

ALARMS

LOADLOCK1


PURGE

REGEN

F1 F2 F3 F4 F5 F6

Confirm gas is properly hooked up with a regulator capable of supplying 60 PSI of gas pressure at a flow rate of 200 l/m.

- From the main screen:
 - Verify Chamber Pressure Control is **ON**, **ON** should be displayed.
 - Chamber Circulation should be **OFF**, **OFF** should be displayed.



LC TECHNOLOGY SOLUTIONS EL103

CHAMBER PRESSURE: +0000.0 mBAR

CHAMBER OXYGEN LEVEL: +0000.0 PPM

CHAMBER MOISTURE LEVEL: +0000.0 PPM

OFF

CHAMBER LIGHTS

OFF

CHAMBER PRESSURE CONTROL

CHAMBER CIRCULATION

CHAMBER SETTINGS

ALARMS

LOADLOCK1

PURGE

REGEN

F1 F2 F3 F4 F5 F6


- From the main menu:
 - Press the Purge button.



LC TECHNOLOGY SOLUTIONS INC.

2C Fanaras Drive, Salisbury, MA 01952 • Phone: 978-255-1620 • Fax: 978-428-0222

E-mail: info@lctechinc.com • Web: www.lctechinc.com



CHAMBER PURGE CONTROL:

ON CHAMBER PURGE TIME SP: 000 MIN


CHAMBER PURGE CHAMBER PURGE TIME REMAINING: 000 : 00

< RETURN

F1 F2 F3 F4 F5 F6

3. Press the Chamber Purge button **ON** will be displayed and the system will begin to purge automatically.

4. The Automatic Purge Function is controlled by time. The purge time is preset at the factory for twenty (20) minutes.



CHAMBER PURGE CONTROL:

OFF CHAMBER PURGE TIME SP: 000 MIN

CHAMBER PURGE CHAMBER PURGE TIME REMAINING: 000 : 00

< RETURN

F1 F2 F3 F4 F5 F6

5. From the Chamber Purge Control screen check the purge time.

To change the length of time follow these steps:

- Press the white box to the right of Chamber Purge Time SP: MIN
- Enter desired set point
- Press Enter
- Press Return

NOTE: Do not exceed (45) minutes if purging with cylinders.

For dewars and house gas supplies refer to Section 12 for purge times.



LC TECHNOLOGY SOLUTIONS INC.

2C Fanaras Drive, Salisbury, MA 01952 • Phone: 978-255-1620 • Fax: 978-428-0222

E-mail: info@lctechinc.com • Web: www.lctechinc.com

CHAMBER PURGE CONTROL:

OFF

CHAMBER PURGE

CHAMBER PURGE TIME SP: 000 MIN

CHAMBER PURGE TIME REMAINING: 000 : 00

< RETURN

F1

F2

F3

F4

F5

F6

6. Chamber Purge Time
Remaining reflects the amount of time remaining in the purge cycle.


NOTE: Once the purge cycle time ends the system will stop purging.



LC TECHNOLOGY SOLUTIONS INC.

2C Fanaras Drive, Salisbury, MA 01952 • Phone: 978-255-1620 • Fax: 978-428-0222
E-mail: info@lctechinc.com • Web: www.lctechinc.com

13.4 Automatic Purge Function Connected to Oxygen Level




This option allows the user to set the automatic purge function to the alarm levels of the oxygen and moisture inside the glovebox. Most customers use this to automatically purge the glovebox if the oxygen level exceeds the alarm.

1. Set the desired oxygen alarm level on the setting screen.


NOTE: This is typically set to 10 ppm.

NOTE: To deactivate this function set the alarm level to 1000 ppm.



When the oxygen level reaches the alarm set point the system turns off the circulation and starts to purge the system.

It will continue to purge until the oxygen level is below the alarm set point.




When the alarm is triggered the automatic purge will activate and run for the amount of time set on the Chamber Purge Control screen. It will continue to purge even after the alarm has cleared for the time set in this function screen.




LC TECHNOLOGY SOLUTIONS INC.

2C Fanaras Drive, Salisbury, MA 01952 • Phone: 978-255-1620 • Fax: 978-428-0222
E-mail: info@lctechinc.com • Web: www.lctechinc.com



When the purge time has elapsed the system will automatically turn circulation on and purge off.

If the alarm triggers again, it will repeat the above steps until the system is able to maintain an oxygen level below the alarm set point.



If the system was purchased with an alarm light and/or buzzer, the alarm/buzzer will trigger when the system is in an alarm state.

The buzzer has an adjustable volume control lever setting on the front of the light tower.




LC TECHNOLOGY SOLUTIONS INC.

2C Fanaras Drive, Salisbury, MA 01952 • Phone: 978-255-1620 • Fax: 978-428-0222
E-mail: info@lctechinc.com • Web: www.lctechinc.com


13.5 Circulation

Before turning on circulation it is very important to make sure the glovebox system has been purged properly.



The screenshot shows the LC Technology Solutions EL103 control panel. At the top, it displays 'CHAMBER PRESSURE: +0000.0 mBAR', 'CHAMBER OXYGEN LEVEL: +0000.0 PPM', and 'CHAMBER MOISTURE LEVEL: +0000.0 PPM'. Below these are buttons for 'OFF', 'CHAMBER LIGHTS', 'ON', 'CHAMBER PRESSURE CONTROL', 'CHAMBER CIRCULATION', and 'CHAMBER SETTINGS'. The 'CHAMBER CIRCULATION' button is highlighted with a yellow circle. To the right of the main panel are buttons for 'ALARMS', 'LOADLOCK1', 'PURGE', and 'REGEN'. At the bottom are function keys F1 through F6.

1. To start circulation of the glovebox environment through the filter column, Chamber Pressure Control must be **ON**.
2. To begin circulation press Chamber Circulation. **ON** will be displayed.




The screenshot shows the LC Technology Solutions EL103 control panel. At the top, it displays 'CHAMBER PRESSURE: +0000.0 mBAR', 'CHAMBER OXYGEN LEVEL: +0000.0 PPM', and 'CHAMBER MOISTURE LEVEL: +0000.0 PPM'. Below these are buttons for 'OFF', 'CHAMBER LIGHTS', 'OFF', 'CHAMBER PRESSURE CONTROL', 'CHAMBER CIRCULATION', and 'CHAMBER SETTINGS'. The 'CHAMBER CIRCULATION' button is highlighted with a yellow circle. To the right of the main panel are buttons for 'ALARMS', 'LOADLOCK1', 'PURGE', and 'REGEN'. At the bottom are function keys F1 through F6.

3. To turn off circulation, press Chamber Circulation. **OFF** will be displayed.




Section 14: Antechamber Operation


14.1 Bringing Items into the Glovebox




1. Make sure antechamber is filled and the vacuum gauge reads zero.
2. Open outside antechamber door.



3. Load green bin or sliding tray with desired material. If using a sliding tray place items on the tray nearest to the inside antechamber door.



4. Place green bin inside antechamber; bin should be loaded into chamber short side first.




5. Push the bin all the way into antechamber until it touches inside antechamber door.

If using a system with a sliding tray slide tray all the way back into chamber.



LC TECHNOLOGY SOLUTIONS INC.


2C Fanaras Drive, Salisbury, MA 01952 • Phone: 978-255-1620 • Fax: 978-428-0222
E-mail: info@lctechinc.com • Web: www.lctechinc.com



6. Close outside antechamber door.

CAUTION: Do not over tighten.

Manual Evacuation



Evacuation Hand Valve


7. For glovebox models with manual evacuation valves turn the evacuation hand valve, located on top of antechamber, to the left.

Antechamber will begin to evacuate.

Continue to evacuate until the vacuum gauge reads -30.

OR

From Main Screen press Loadlock. Loadlock Control screen will now be displayed.




Press the Evacuation button on glovebox systems with touch screen controls to evacuate the antechamber.



LC TECHNOLOGY SOLUTIONS INC.


2C Fanaras Drive, Salisbury, MA 01952 • Phone: 978-255-1620 • Fax: 978-428-0222

E-mail: info@lctechinc.com • Web: www.lctechinc.com



8. Close evacuation valve by turning the evacuation hand valve to the right.

OR



Push the EVAC button to stop evacuation.



LC TECHNOLOGY SOLUTIONS INC.

2C Fanaras Drive, Salisbury, MA 01952 • Phone: 978-255-1620 • Fax: 978-428-0222

E-mail: info@lctechinc.com • Web: www.lctechinc.com

Swagelok Refill Valve




9. For glovebox models with manual refill, refill antechamber using Swagelok refill valve.

Continue to refill until the vacuum gauge reads 0.

OR

Press REFILL button on glovebox systems with automatic controls to begin process. Push the REFILL button to stop.



10. Repeat the above cycle two (2) more times for a total of three (3) evacuation/refill cycles.

NOTE: Before opening the inside door make sure vacuum gauge reads zero (0) & evacuation and refill valves are closed/turned off.



LC TECHNOLOGY SOLUTIONS INC.

2C Fanaras Drive, Salisbury, MA 01952 • Phone: 978-255-1620 • Fax: 978-428-0222
E-mail: info@lctechinc.com • Web: www.lctechinc.com



11. Open inside door, remove green bin, and close antechamber door.

NOTE: When opening the inside door, spin door handle completely until door makes contact with door arm. This is important so door does not make contact with side wall of glovebox and damage the sealing surface.



14.2 Removing Items from Glovebox

Determine whether the antechamber has room air or inert gas in it.

If inert gas follow these steps:

1. Open inside antechamber door.
2. Load green bin/tray into antechamber.
3. Close inside antechamber door.
4. Open outside antechamber door.

If room air follow these steps:

1. Evacuate and refill the antechamber three (3) times.
2. Refer to Section 14.1, Steps 6-9 for the evacuation/refill process.


NOTE: This will ensure the antechamber has inert gas in it.

1. Open inside antechamber door.
2. Load green bin/tray into antechamber.
3. Close inside antechamber door.
4. Open outside antechamber door.



14.3 Automatic Antechamber Control / Loadlock

The automatic antechamber control / loadlock function is used to automatically cycle the antechamber between the evacuation function and the refill function. The cycle is controlled by the number of cycles selected and the vacuum level selected.



LC TECHNOLOGY SOLUTIONS EL103


CHAMBER PRESSURE: +0000.0 mBAR
CHAMBER OXYGEN LEVEL: +0000.0 PPM
CHAMBER MOISTURE LEVEL: +0000.0 PPM

ALARMS
LOADLOCK1
PURGE
REGEN
CHAMBER LIGHTS
CHAMBER PRESSURE CONTROL
CHAMBER CIRCULATION
CHAMBER SETTINGS

F1 F2 F3 F4 F5 F6

1. Press the Loadlock button from the main screen.

Loadlock Control screen will be displayed.



LOADLOCK1 CONTROL:

REFILL
OFF
EVACUATE
< RETURN

EVACUATION PRESSURE SP: +0000.00 mBAR
REFILL TIME SP: +00000 MIN
LAST REFILL TIME SP: +00000 MIN
EVACUATION CYCLES SP: +00000
EVACUATION PRESSURE: +0000.00 mBAR
OPERATING MODE: MANUAL
AUTO CYCLE: OFF
AUTO CYCLE STATUS: OFF
REFILL TIME REMAINING: 000 : 00
CYCLES REMAINING: 000

F1 F2 F3 F4 F5 F6

Once all of the settings below have been set, press the Evacuate button to begin the automatic cycle.



LC TECHNOLOGY SOLUTIONS INC.

2C Fanaras Drive, Salisbury, MA 01952 • Phone: 978-255-1620 • Fax: 978-428-0222

E-mail: info@lctechinc.com • Web: www.lctechinc.com

LOADLOCK1 CONTROL:

OFF EVACUATION PRESSURE SP: +0000.00 mBAR

REFILL REFILL TIME SP: +00000 MIN

LAST REFILL TIME SP: +00000 MIN

OFF EVACUATION CYCLES SP: +00000

EVACUATE EVACUATION PRESSURE: +0000.00 mBAR

OPERATING MODE: MANUAL

AUTO CYCLE: OFF

AUTO CYCLE STATUS: OFF

REFILL TIME REMAINING: 000 : 00

< RETURN CYCLES REMAINING: 000

F1 F2 F3 F4 F5 F6

2. The automatic antechamber control / loadlock settings can be found on this screen.
3. To change these set points, press the white box to the right of desired setting. Type in new setting and press Enter.

LOADLOCK1 CONTROL:

OFF EVACUATION PRESSURE SP: +0000.00 mBAR

REFILL REFILL TIME SP: +00000 MIN

LAST REFILL TIME SP: +00000 MIN

OFF EVACUATION CYCLES SP: +00000

EVACUATE EVACUATION PRESSURE: +0000.00 mBAR

OPERATING MODE: MANUAL

AUTO CYCLE: OFF

AUTO CYCLE STATUS: OFF

REFILL TIME REMAINING: 000 : 00

< RETURN CYCLES REMAINING: 000

F1 F2 F3 F4 F5 F6

There are four settings that can be changed in the Loadlock1 Control screen.

1. Evacuation Pressure SP - This is the evacuation level that will be reached before refilling.

NOTE: 0.5 mBAR is the recommended Evacuation Pressure SP.



LC TECHNOLOGY SOLUTIONS INC.

2C Fanaras Drive, Salisbury, MA 01952 • Phone: 978-255-1620 • Fax: 978-428-0222
E-mail: info@lctechinc.com • Web: www.lctechinc.com

LOADLOCK1 CONTROL:

OFF

REFILL

OFF

EVACUATE

< RETURN

EVACUATION PRESSURE SP: +0000.00 mBAR

REFILL TIME SP: +00000 MIN

LAST REFILL TIME SP: +00000 MIN

EVACUATION CYCLES SP: +00000

EVACUATION PRESSURE: +0000.00 mBAR

OPERATING MODE: MANUAL

AUTO CYCLE: OFF

AUTO CYCLE STATUS: OFF

REFILL TIME REMAINING: 000 : 00

CYCLES REMAINING: 000

F1

F2

F3

F4

F5

F6

2. Refill Time SP - This is the length of time the chamber will refill before another evacuation cycle is started.

NOTE: +1 Minute is the recommended Refill Time SP.

LOADLOCK1 CONTROL:

OFF

REFILL

OFF

EVACUATE

< RETURN

EVACUATION PRESSURE SP: +0000.00 mBAR

REFILL TIME SP: +00000 MIN

LAST REFILL TIME SP: +00000 MIN

EVACUATION CYCLES SP: +00000

EVACUATION PRESSURE: +0000.00 mBAR

OPERATING MODE: MANUAL

AUTO CYCLE: OFF

AUTO CYCLE STATUS: OFF

REFILL TIME REMAINING: 000 : 00

CYCLES REMAINING: 000

F1

F2

F3

F4

F5

F6

3. Last Refill Time SP - This set point is the final refill cycle and should be longer than the second refill time.

NOTE: +2 Minutes is the recommended Last Refill Time SP.

LOADLOCK1 CONTROL:

OFF

REFILL

OFF

EVACUATE

< RETURN

EVACUATION PRESSURE SP: +0000.00 mBAR

REFILL TIME SP: +00000 MIN

LAST REFILL TIME SP: +00000 MIN

EVACUATION CYCLES SP: +00000

EVACUATION PRESSURE: +0000.00 mBAR

OPERATING MODE: MANUAL

AUTO CYCLE: OFF

AUTO CYCLE STATUS: OFF

REFILL TIME REMAINING: 000 : 00

CYCLES REMAINING: 000

F1

F2

F3

F4

F5

F6

4. Evacuation Cycles SP - This allows the user to set the number of EVAC/Refill cycles the chamber will perform.

NOTE: The recommended number of cycles is three (3).



LC TECHNOLOGY SOLUTIONS INC.

2C Fanaras Drive, Salisbury, MA 01952 • Phone: 978-255-1620 • Fax: 978-428-0222

E-mail: info@lctechinc.com • Web: www.lctechinc.com

LOADLOCK1 CONTROL:

OFF

REFILL

OFF

EVACUATE

< RETURN

EVACUATION PRESSURE SP: +0000.00 mBAR

REFILL TIME SP: +00000 MIN

LAST REFILL TIME SP: 00000 MIN

EVACUATION CYCLES SP: +00000

EVACUATION PRESSURE: -0000.00 mBAR

OPERATING MODE: MANUAL

AUTO CYCLE: OFF

AUTO CYCLE STATUS: OFF

REFILL TIME REMAINING: 000 : 00

CYCLES REMAINING: 000

F1

F2

F3

F4

F5

F6

- To begin automatic antechamber control change the Evacuation Cycles SP: to 3 by pressing the white button to the right.
- Enter the set point and press Return.




The Operating Mode will change from Manual to Auto.

- Next press the Evacuate button to begin the automatic cycle.

During the automatic cycle the following information can be found on the Loadlock1 Control screen.

- a. Auto Cycle: ON – This is the indicator that the auto cycle is running.
- b. Auto Cycle Status: Evac/Refill
– This status indicates if the chamber is evacuating or refilling.
- c. Refill Time Remaining – This status indicates the amount of time remaining in the refill process.
- d. Cycles Remaining – This status indicates the number of cycles remaining. The number of cycles will count down from 3, 2, 1, 0.



The vacuum level is shown on the gauge mounted on top of the antechamber and on the PLC screen.

Once the cycle is complete the evacuation and refill cycles will turn off.



LC TECHNOLOGY SOLUTIONS INC.

2C Fanaras Drive, Salisbury, MA 01952 • Phone: 978-255-1620 • Fax: 978-428-0222

E-mail: info@lctechinc.com • Web: www.lctechinc.com

LOADLOCK1 CONTROL:

OFF EVACUATION PRESSURE SP: +0000.00 mBAR

REFILL REFILL TIME SP: +00000 MIN

LAST REFILL TIME SP: +00000 MIN

OFF EVACUATION CYCLES SP: +00000

EVACUATE EVACUATION PRESSURE: +0000.00 mBAR

OPERATING MODE: **MANUAL**

AUTO CYCLE: OFF

AUTO CYCLE STATUS: OFF

REFILL TIME REMAINING: 000 : 00

< RETURN CYCLES REMAINING: 000

F1 F2 F3 F4 F5 F6

NOTE: System can also be operated in manual mode by setting Evacuation Cycles SP to zero (0).

Operating Mode will now display Manual. Refer to Section 14.1 steps 6-9 for operational instructions.



LC TECHNOLOGY SOLUTIONS INC.


2C Fanaras Drive, Salisbury, MA 01952 • Phone: 978-255-1620 • Fax: 978-428-0222
E-mail: info@lctechinc.com • Web: www.lctechinc.com

14.4 Manual Purge for Large Antechamber


1. Open outer door and load large antechamber.



2. Close door.



3. Open vent valve.



4. Open gas valve.

NOTE: Connect inert gas to gas valve. (5 psi maximum for purging.)
NOTE: Automatic antechamber control cycles needs to be set to Zero so that the door locks are disabled during this process. See page 44.

5. Purge for 10 minutes.

6. Close gas valve and vent valve.

7. Open inside large antechamber door.



LC TECHNOLOGY SOLUTIONS INC.

2C Fanaras Drive, Salisbury, MA 01952 • Phone: 978-255-1620 • Fax: 978-428-0222
E-mail: info@lctechinc.com • Web: www.lctechinc.com

14.5 Mini Antechamber Operation




1. Open outside mini antechamber door by pulling up on red lever.

2. Remove outside mini antechamber door.




3. Load green bin with parts.
4. Slide green bin all the way into antechamber until it comes into contact with inside antechamber door.

5. Replace outside antechamber door and push red lever down into the locked position.



6. Turn hand valve to evacuation position.
7. Continue evacuating until gauge reaches approximately -30.




8. Turn hand valve to refill position.
9. Continue to refill until the vacuum gauge reads 0.



LC TECHNOLOGY SOLUTIONS INC.

2C Fanaras Drive, Salisbury, MA 01952 • Phone: 978-255-1620 • Fax: 978-428-0222
E-mail: info@lctechinc.com • Web: www.lctechinc.com

10. Repeat the above cycle two (2) more times for a total of three (3) evacuation/refill cycles.



11. Return hand valve to up (closed) position.

12. Remove inside antechamber door.

13. Remove green bin and return to chamber when complete.

14. Replace inside antechamber door and push red lever into the locked position.

15. Refer to Section 14.2 for instructions for Removing Items from Glovebox.




LC TECHNOLOGY SOLUTIONS INC.

2C Fanaras Drive, Salisbury, MA 01952 • Phone: 978-255-1620 • Fax: 978-428-0222
E-mail: info@lctechinc.com • Web: www.lctechinc.com


14.6 Mini Antechamber with Automatic Control for Door Lock




1. To operate mini antechamber press LoadLock2 button.



2. Push Outer Door button.



3. Open door and load chamber.




4. Close door.




LC TECHNOLOGY SOLUTIONS INC.

2C Fanaras Drive, Salisbury, MA 01952 • Phone: 978-255-1620 • Fax: 978-428-0222

E-mail: info@lctechinc.com • Web: www.lctechinc.com




5. Push Outer Door button.



6. Push Evacuate button. This will allow it to run for 3 cycles.

NOTE: For setting parameters refer to Section: 14.3 Automatic Antechamber Control / Loadlock.



NOTE: You can push the Refill button any time to abort a cycle.




LC TECHNOLOGY SOLUTIONS INC.


2C Fanaras Drive, Salisbury, MA 01952 • Phone: 978-255-1620 • Fax: 978-428-0222
E-mail: info@lctechinc.com • Web: www.lctechinc.com

14.7 Manual Purge for Mini Antechamber


1. Open outer door and load mini antechamber.



2. Close door.



3. Open vent valve.



4. Open gas valve.

NOTE: Connect inert gas to gas valve. (5 psi maximum for purging.)
NOTE: Automatic antechamber control cycles needs to be set to Zero so that the door locks are disabled during this process. See page 44.

5. Purge for 5 minutes.

6. Close gas valve and vent valve.


7. Open inside mini antechamber door.



Section 15: Regeneration Mode

Regeneration Mode is for reactivating the filter material. It can only be activated if Chamber Circulation is **OFF** and Chamber Pressure is **ON**.


From the main screen:



Prior to running a regeneration if Oxygen levels are greater than 250 PPM purge the glovebox for 20 minutes.

1. Turn Circulation **OFF**
2. Push REGEN Button.

The screen will show you the following message:



3. REGEN MODE will now display.


Caution: The regeneration cycle times are critical to the system running properly. Do not change these settings!

It is recommended that you contact LC Technology at (978) 255-1620 before making any changes.




LC TECHNOLOGY SOLUTIONS INC.

2C Fanaras Drive, Salisbury, MA 01952 • Phone: 978-255-1620 • Fax: 978-428-0222
E-mail: info@lctechinc.com • Web: www.lctechinc.com



1. To start a regeneration press the Gas Purifier Regeneration button. **ON** will be displayed.
2. To change the set points for heat, purge, evac, or cool press the appropriate white box to the right of the Time SP. Enter the number of minutes for the cycle and press Return.

The screen will show you the following message:




3. The following Message will appear, make sure regeneration gas is connected to the system.
4. Once you have connected your regeneration gas (**NOTE: the regen gas should be set at 15 PSI, with a flow rate of 15L per minute**) you need to confirm this message by pressing the Acknowledge Button.
5. Regeneration will start automatically and the heat time will start to run.



LC TECHNOLOGY SOLUTIONS INC.

2C Fanaras Drive, Salisbury, MA 01952 • Phone: 978-255-1620 • Fax: 978-428-0222

E-mail: info@lctechinc.com • Web: www.lctechinc.com



During the regeneration cycle the following information can be found on the Regen screen.

- A. Time SP – This indicate the number of minutes for each process, heat, purge, evac, and cool.
- B. Remaining – This is the number of minutes remaining in each cycle of heating, purging, evacuating and cooling.

The regeneration process takes (13) hours and has (4) phases. Each phase will be displayed as the program progresses from one phase to the other.

- The first phase is heating which lasts three (3) hours depending on the system.
 - **NOTE:** You may smell a slight odor at this time. This is normal.
- The second phase is purging which lasts (3) hours.
 - **NOTE:** Confirm the regeneration gas flow rate at this time at 15L/min. Adjust if needed at the regeneration gas flow meter.




- The third phase is evacuation which lasts (3) hours.
- The fourth phase is cooling which lasts (4) hours.



LC TECHNOLOGY SOLUTIONS INC.

2C Fanaras Drive, Salisbury, MA 01952 • Phone: 978-255-1620 • Fax: 978-428-0222
E-mail: info@lctechinc.com • Web: www.lctechinc.com




After 13 hours the Regeneration Cycle will be complete:

6. Press Gas Purifier Regeneration Button to turn off Regeneration Function, **OFF** will be displayed.
7. Return to main screen and restart circulation by Pressing Circulation **ON** should be displayed.

To interrupt the regeneration process, press the Gas Purifier Regeneration button at any time. Caution should be used at this time.

- If the system has been heating for over (1) hour, the system should be allowed to cool for at least (2) hours before the regeneration is restarted.
- If the system has past the heating phase, the system should cool for (6) hours before restarting the regeneration.

Please call LC Technology at (978) 255-1620 if you have any questions about interrupting the regeneration or restarting the regeneration process.




This screen shows that the regeneration has been interrupted due to the system being shut down during the regeneration cycle.

Pressing the Continue button allows the cycle to continue.

Pressing the Abort button exits the interrupted cycle and the regeneration mode.



15.1 Testing the Regeneration Process



IMPORTANT:

To test the regeneration process follow the steps below:

1. After the regeneration process has completed empty the drain portion of the common vent line by removing the capping nut that is located at the bottom of the common vent pipe.
2. Measure the amount of water collected.

NOTE: A normal regeneration yields 50 milliliters of water. If the water is less it could mean something is wrong.

NOTE: It is recommended to change your vacuum pump oil after a regeneration.



LC TECHNOLOGY SOLUTIONS INC.

2C Fanaras Drive, Salisbury, MA 01952 • Phone: 978-255-1620 • Fax: 978-428-0222
E-mail: info@lctechinc.com • Web: www.lctechinc.com

15.2 Regeneration of a Dual Filter Column System

CHAMBER SETTINGS:

CHAMBER PRESSURE CONTROL HIGH SP:	+0000.0	mBAR
CHAMBER PRESSURE CONTROL LOW SP:	+0000.0	mBAR
CHAMBER GAS SOLENOID OFF DELAY SP:	+00000	x.1SEC
CHAMBER VACUUM SOLENOID OFF DELAY SP:	+00000	x.1SEC
CHAMBER PRESSURE HIGH ALARM SP:	+0000.0	mBAR
CHAMBER PRESSURE LOW ALARM SP:	+0000.0	mBAR
CHAMBER OXYGEN LEVEL HIGH ALARM SP:	+0000.0	PPM
CHAMBER MOISTURE LEVEL HIGH ALARM SP:	+0000.0	PPM

< RETURN CHAMBER SELECT CHAMBER 1 & 2

FILTER SELECT FILTER 1

F1 F2 F3 F4 F5 F6

If your system was configured with dual filter columns you can circulate your system and regenerate the other filter column at the same time. The system automatically regenerates the filter column that is not being used for circulation mode. In the chamber settings menu, filter 1 is selected for circulation therefore filter 2 will be regenerated. Please follow the instruction in Section 15 to regenerate your filter column.



Section 16: Chamber Settings

16.1 Set Points

All user and system set points are located in Chamber Settings. These set points have been preset at the factory and care should be taken when altering the set points.

To enter Chamber Settings:

1. Press Chamber Settings button.
2. Select Set Point you wish to change.
3. Enter new Set Point.

In Set Points you can change the following settings:

1. Chamber Pressure Control High/Low SP
2. Chamber Gas / Vac Solenoid Delay SP
3. Chamber Pressure High/Low Alarm SP
4. Chamber Oxygen Level High Alarm SP
5. Chamber Moisture Level High Alarm SP

To alter the set points press desired Set Point and enter new Set Point.

Chamber Pressure Control High & Low Set Points

CHAMBER SETTINGS:

CHAMBER PRESSURE CONTROL HIGH SP:	+0000.0	mBAR
CHAMBER PRESSURE CONTROL LOW SP:	+0000.0	mBAR
CHAMBER GAS SOLENOID OFF DELAY SP:	+00000	x.1SEC
CHAMBER VACUUM SOLENOID OFF DELAY SP:	+00000	x.1SEC
CHAMBER PRESSURE HIGH ALARM SP:	+0000.0	mBAR
CHAMBER PRESSURE LOW ALARM SP:	+0000.0	mBAR
CHAMBER OXYGEN LEVEL HIGH ALARM SP:	+0000.0	PPM
CHAMBER MOISTURE LEVEL HIGH ALARM SP:	+0000.0	PPM

< RETURN CHAMBER SELECT CHAMBER 1 & 2

F1 F2 F3 F4 F5 F6

This screen allows the user to select the box pressure ranges that are needed for your application. Typical settings are 0.5 mBAR for low limit and 2.5 mBAR for high limit.

The system will maintain the pressure in the glovebox between these two values. The foot pedals will only adjust the pressure in the glovebox between these ranges.



LC TECHNOLOGY SOLUTIONS INC.

2C Fanaras Drive, Salisbury, MA 01952 • Phone: 978-255-1620 • Fax: 978-428-0222

E-mail: info@lctechinc.com • Web: www.lctechinc.com

Chamber Gas & Vac Solenoid Off Delay Set Points

CHAMBER SETTINGS:

CHAMBER PRESSURE CONTROL HIGH SP:	+0000.0	mBAR
CHAMBER PRESSURE CONTROL LOW SP:	+0000.0	mBAR
CHAMBER GAS SOLENOID OFF DELAY SP:	+00000	x.1SEC
CHAMBER VACUUM SOLENOID OFF DELAY SP:	+00000	x.1SEC
CHAMBER PRESSURE HIGH ALARM SP:	+0000.0	mBAR
CHAMBER PRESSURE LOW ALARM SP:	+0000.0	mBAR
CHAMBER OXYGEN LEVEL HIGH ALARM SP:	+0000.0	PPM
CHAMBER MOISTURE LEVEL HIGH ALARM SP:	+0000.0	PPM

< RETURN CHAMBER SELECT CHAMBER 1 & 2

F1 F2 F3 F4 F5 F6

This set point allows the user the ability to set a delay on how long the gas and/or vac valves stay open when triggered.

The typical setting is 5 for the Chamber Gas Solenoid Off Delay SP and the Chamber Vacuum Solenoid Off Delay SP.

Chamber Pressure High & Low Alarm Set Points

CHAMBER SETTINGS:

CHAMBER PRESSURE CONTROL HIGH SP:	+0000.0	mBAR
CHAMBER PRESSURE CONTROL LOW SP:	+0000.0	mBAR
CHAMBER GAS SOLENOID OFF DELAY SP:	+00000	x.1SEC
CHAMBER VACUUM SOLENOID OFF DELAY SP:	+00000	x.1SEC
CHAMBER PRESSURE HIGH ALARM SP:	+0000.0	mBAR
CHAMBER PRESSURE LOW ALARM SP:	+0000.0	mBAR
CHAMBER OXYGEN LEVEL HIGH ALARM SP:	+0000.0	PPM
CHAMBER MOISTURE LEVEL HIGH ALARM SP:	+0000.0	PPM

< RETURN CHAMBER SELECT CHAMBER 1 & 2

F1 F2 F3 F4 F5 F6

Chamber Pressure High Alarm Set Point: This allows the user to input the high pressure alarm set point. The factory default is +15.0 mBAR.

Chamber Pressure Low Alarm Set Point: This allows the user to input the low pressure alarm set point. The factory default is -15.0 mBAR.



LC TECHNOLOGY SOLUTIONS INC.

2C Fanaras Drive, Salisbury, MA 01952 • Phone: 978-255-1620 • Fax: 978-428-0222

E-mail: info@lctechinc.com • Web: www.lctechinc.com

Chamber Oxygen & Moisture Level High Alarm Set Points

This allows the user to input the oxygen and moisture alarm set points. The factory default is 50 ppm.

CHAMBER SETTINGS:

CHAMBER PRESSURE CONTROL HIGH SP:	+0000.0	mBAR
CHAMBER PRESSURE CONTROL LOW SP:	+0000.0	mBAR
CHAMBER GAS SOLENOID OFF DELAY SP:	+00000	x.1SEC
CHAMBER VACUUM SOLENOID OFF DELAY SP:	+00000	x.1SEC
CHAMBER PRESSURE HIGH ALARM SP:	+0000.0	mBAR
CHAMBER PRESSURE LOW ALARM SP:	+0000.0	mBAR
CHAMBER OXYGEN LEVEL HIGH ALARM SP:	+0000.0	PPM
CHAMBER MOISTURE LEVEL HIGH ALARM SP:	+0000.0	PPM

< RETURN

CHAMBER SELECT CHAMBER 1 & 2

F1 F2 F3 F4 F5 F6



LC TECHNOLOGY SOLUTIONS INC.

2C Fanaras Drive, Salisbury, MA 01952 • Phone: 978-255-1620 • Fax: 978-428-0222

E-mail: info@lctechinc.com • Web: www.lctechinc.com

16.2 Enable / Disable

Chamber Selection

CHAMBER SETTINGS:

CHAMBER PRESSURE CONTROL HIGH SP:	+0000.0	mBAR
CHAMBER PRESSURE CONTROL LOW SP:	+0000.0	mBAR
CHAMBER GAS SOLENOID OFF DELAY SP:	+00000	x.1SEC
CHAMBER VACUUM SOLENOID OFF DELAY SP:	+00000	x.1SEC
CHAMBER PRESSURE HIGH ALARM SP:	+0000.0	mBAR
CHAMBER PRESSURE LOW ALARM SP:	+0000.0	mBAR
CHAMBER OXYGEN LEVEL HIGH ALARM SP:	+0000.0	PPM
CHAMBER MOISTURE LEVEL HIGH ALARM SP:	+0000.0	PPM

< RETURN

CHAMBER SELECT

CHAMBER 1 & 2

F1

F2

F3

F4

F5

F6

1. Customers with 2 separate boxes that are run off of one purifier cart are capable of disabling or enabling one or both boxes by using bypass isolation valve.
2. Press Chamber Select to select desired chamber.


Note: You may only have one box.



LC TECHNOLOGY SOLUTIONS INC.

2C Fanaras Drive, Salisbury, MA 01952 • Phone: 978-255-1620 • Fax: 978-428-0222
E-mail: info@lctechinc.com • Web: www.lctechinc.com

16.3 Dual Filter Columns



CHAMBER SETTINGS:

CHAMBER PRESSURE CONTROL HIGH SP:	+0000.0	mBAR
CHAMBER PRESSURE CONTROL LOW SP:	+0000.0	mBAR
CHAMBER GAS SOLENOID OFF DELAY SP:	+00000	x.1SEC
CHAMBER VACUUM SOLENOID OFF DELAY SP:	+00000	x.1SEC
CHAMBER PRESSURE HIGH ALARM SP:	+0000.0	mBAR
CHAMBER PRESSURE LOW ALARM SP:	+0000.0	mBAR
CHAMBER OXYGEN LEVEL HIGH ALARM SP:	+0000.0	PPM
CHAMBER MOISTURE LEVEL HIGH ALARM SP:	+0000.0	PPM

< RETURN CHAMBER SELECT CHAMBER 1 & 2

FILTER SELECT FILTER 1

F1 F2 F3 F4 F5 F6

1. For systems supplied with 2 filter columns, you will need to select the active filter column. The activate filter column is the column that you wish to use to circulate the glovebox environment. The other column will be in standby mode.
2. Press Filter Select to select desired filter. (Filter 1 or Filter 2).

Do not switch filter column while running a regeneration.




LC TECHNOLOGY SOLUTIONS INC.

2C Fanaras Drive, Salisbury, MA 01952 • Phone: 978-255-1620 • Fax: 978-428-0222
E-mail: info@lctechinc.com • Web: www.lctechinc.com

Section 17: Analyzers

The display will read as follows:



The main screen provides the current PPM level of oxygen and moisture inside the chamber.




Section 18: Solvent Removal Systems, Operation and Maintenance


18.1 Manual Solvent Removal System Operation

The solvent removal system is for the removal of solvent vapors from the glovebox environment. The system has two (2) operation modes: active mode and bypass mode.

NOTE: There are three (3) main valves on the top of the solvent removal system. The valves are labeled 1, 2 and 3. It is very important these valves are not all closed at the same time. An open path from the gas purification system to the glovebox must be maintained at all times, otherwise serious damage may occur to the system.



Valve 1




Valve 2 and 3

Active Mode:

This is the mode the system will normally be in. This will allow the glovebox environment to circulate through the solvent removal system so it can trap solvent vapor.

To put the system into Active Mode, valve 1 should already be open. Open valves 2 and 3 and close valve 1.



Valve 1



Valve 2 and 3

Bypass Mode:

This mode is only used for maintaining the solvent removal system. The system should only be put in Bypass Mode so that the activated carbon in the solvent removal system can be changed.

Open valve 1 and close valves 2 and 3.

To put the system back into Active Mode, open valves 2 and 3, then close valve 1.




18.2 Manual Solvent Removal System Maintenance

Approximately every 3-6 months the activated carbon in the solvent removal system will have to be changed.

Follow the steps below to change the activated carbon:

1. Put the solvent removal system in Bypass Mode.




2. Remove KF40 clamp and cover on the fill/empty port on top of solvent removal system.



3. Using a shop vac and solvent extraction tool suck all of the used activated carbon out of the solvent trap.

4. Using a funnel refill the solvent removal system with fresh activated carbon. It will hold 10 lbs of material. Do not fill with more than this amount.



5. Replace KF40 cover and clamp on top of solvent filter.



6. Using Valve 4 mounted on the back of the system evacuate the solvent removal system for approximately 24 hrs by turning Valve 4 to the left.



LC TECHNOLOGY SOLUTIONS INC.

2C Fanaras Drive, Salisbury, MA 01952 • Phone: 978-255-1620 • Fax: 978-428-0222
E-mail: info@lctechinc.com • Web: www.lctechinc.com

7. Using Valve 4 refill the solvent trap with inert gas by turning Valve 4 to the right. The solvent trap will fill with inert gas from the glovebox.

8. Return Valve 4 to the center position.

9. Put the system back into active mode, as described on page 62.





LC TECHNOLOGY SOLUTIONS INC.

2C Fanaras Drive, Salisbury, MA 01952 • Phone: 978-255-1620 • Fax: 978-428-0222

E-mail: info@lctechinc.com • Web: www.lctechinc.com

18.3 Automatic Solvent Removal System Operation

The solvent removal system is for the removal of solvent vapors from the glovebox environment. The system has two (2) operation modes: active mode and bypass mode.

	<p>To put the Solvent Removal System into Active Mode Press the Solvent Filter button, ACTIVE will be displayed. This will allow the solvent removal system to remove solvents from the glovebox environment.</p>
	<p>To Bypass Solvent Removal filter push Solvent Filter button and it should say BYPASS. The Bypass function is used when you need to reactivate the solvent removal filter material. In by pass mode the system will not be removing solvent from the glovebox environment.</p>






LC TECHNOLOGY SOLUTIONS INC.

2C Fanaras Drive, Salisbury, MA 01952 • Phone: 978-255-1620 • Fax: 978-428-0222
E-mail: info@lctechinc.com • Web: www.lctechinc.com

18.4 Automatic Solvent Reactivation

The reactivation process should be run every three months.

	<ol style="list-style-type: none"> 1. Put system in Bypass by Pressing Solvent Filter button, should say BYPASS
	<ol style="list-style-type: none"> 2. Press Reactivate Solvent Filter button.
	<ol style="list-style-type: none"> 3. Press Solvent Filter Reactivation button. ON should be displayed. This will start the reactivation process.



LC TECHNOLOGY SOLUTIONS INC.

2C Fanaras Drive, Salisbury, MA 01952 • Phone: 978-255-1620 • Fax: 978-428-0222

E-mail: info@lctechinc.com • Web: www.lctechinc.com


The Solvent Reactivation is going to go through 4 phases.

The First is Heating which lasts 3 hours.

The Second Phase is Reactivation which last 2 hours. During the reactivation phase the system will purge inert gas through the solvent removal filter. The process will use approximately 150 cubic feet of gas during the reactivation. The inert gas is supplied by the inert gas connection already on the gas purifier cart. There is no need to connect a separate gas supply for this process.

The Third is Cooling which last 4 hours.


The Fourth is Backfill which last one minute.



After 9 hours the Reactivation Cycle will be complete the system will say **COMPLETE**,

Press the Solvent Filter Reactivation button to turn the process off. **OFF** should be displayed.

Then hit the Return button to go back to the Main Menu



At the Main Menu put the system back in Active Mode by pressing the Solvent Filter button. **ACTIVE** should be displayed.



LC TECHNOLOGY SOLUTIONS INC.


2C Fanaras Drive, Salisbury, MA 01952 • Phone: 978-255-1620 • Fax: 978-428-0222

E-mail: info@lctechinc.com • Web: www.lctechinc.com

To interrupt the reactivation process, press the Solvent Filter Reactivation button at any time. Caution should be used at this time.

- If the system has been heating for over (1) hour, the system should be allowed to cool for at least (2) hours before the reactivation is restarted.
- If the system has past the heating phase, the system should cool for (6) hours before restarting the reactivation.

Please call LC Technology at (978) 255-1620 if you have any questions about interrupting the reactivation or restarting the reactivation process.

	<p>This screen shows that the Reactivation has been interrupted due to the system being shut down during reactivation cycle.</p> <p>Pressing the Continue button allows the cycle to continue.</p> <p>Pressing the Abort button exits the interrupted cycle and the reactivation mode.</p>
--	--



Section 19: Internal Charcoal Trap

The internal charcoal trap is used for trapping small amounts of trace solvent. The internal charcoal trap comes in a set of (12) and should be changed once per month.

NOTE: The filters are labeled 1 – 12, one for each month of the year.

1. Bring new filter into the glovebox.



2. Pull old filter out of the socket and discard.



3. Replace old filter with the new one by pushing it into place in the socket.

NOTE: The new filter goes on the inlet side to the gas purifier.

The correct filter will be labeled with a number.



LC TECHNOLOGY SOLUTIONS INC.

2C Fanaras Drive, Salisbury, MA 01952 • Phone: 978-255-1620 • Fax: 978-428-0222
E-mail: info@lctechinc.com • Web: www.lctechinc.com



NOTE: Repeat every month. Filters are labeled 1-12, one for each month of the year.



Section 20: Freezer Operation and Maintenance

20.1 Freezer Operation



NOTE: Verify freezer is plugged in.

1. Turn freezer on by pressing the switch to the ON position.
2. Turn freezer off by pressing the switch to the OFF position.

NOTE: Freezer is preset to -35°C by the factory.



To adjust the Set Point:

1. Press the second key to the left on the temperature controller.
2. Use the ↑ or ↓ arrows to change set point.
3. Press the second key to the left again to complete.

NOTE: Do not change any other settings on the temperature controller without first contacting LC Technology.

NOTE: The maximum operating temperature is -35°C. The recommended every day operating temperature is -25°C.



LC TECHNOLOGY SOLUTIONS INC.

2C Fanaras Drive, Salisbury, MA 01952 • Phone: 978-255-1620 • Fax: 978-428-0222
E-mail: info@lctechinc.com • Web: www.lctechinc.com

20.2 Freezer Maintenance


The only user serviceable maintenance item for the freezer is cleaning the cooling fins on the compressor. This will keep dust and other particles from building up on the unit.

Freezer maintenance, other than cleaning the cooling fins, should be performed by a qualified refrigerator repair person.




Section 21: Alarm Messages

The alarm levels are set in the Chamber Settings for both the oxygen and moisture analyzers. If the reading on the PLC is above the alarm set point it will display the following messages.



The screenshot shows a green-bordered screen titled "SYSTEM ALARM SCREEN". Below the title, the text "System Alarm High Oxygen and Moisture" is displayed. At the bottom left of the screen is a button labeled "< RETURN". Below the screen are six buttons labeled F1, F2, F3, F4, F5, and F6.

This alarm displays when oxygen and moisture levels are above the alarm set point.



The screenshot shows a green-bordered screen titled "SYSTEM ALARM SCREEN". Below the title, the text "System Alarm High Moisture" is displayed. At the bottom left of the screen is a button labeled "< RETURN". Below the screen are six buttons labeled F1, F2, F3, F4, F5, and F6.


The moisture level high alarm displays when the moisture level in PPM is above the alarm set point.



LC TECHNOLOGY SOLUTIONS INC.

2C Fanaras Drive, Salisbury, MA 01952 • Phone: 978-255-1620 • Fax: 978-428-0222

E-mail: info@lctechinc.com • Web: www.lctechinc.com



The oxygen level high alarm displays when the oxygen level in PPM is above the alarm set point.

Once the level of H₂O and O₂ is below the alarm level the messages will automatically clear.




LC TECHNOLOGY SOLUTIONS INC.

2C Fanaras Drive, Salisbury, MA 01952 • Phone: 978-255-1620 • Fax: 978-428-0222

E-mail: info@lctechinc.com • Web: www.lctechinc.com

Inlet/Outlet Valves Not Open




This alarm means that the electro-pneumatic valves on top of the filter column are not opening properly.

This is almost always caused by low gas pressure or an empty gas cylinder.

Check the gas supply and make sure the system has at least 60 psi going to it.

Inlet/Outlet Valves Not Closed



This means the valves on top of the filter column are not closed and you will be unable to perform the regeneration.

Call LC Technology for more information.




LC TECHNOLOGY SOLUTIONS INC.

2C Fanaras Drive, Salisbury, MA 01952 • Phone: 978-255-1620 • Fax: 978-428-0222

E-mail: info@lctechinc.com • Web: www.lctechinc.com

Over or Under Pressure Alarm



The Low Pressure Alarm means the system is under-pressurized. Once the glovebox goes above the low pressure safety set point it will shut off.

The High Pressure Alarm means the system is over-pressurized. Once the glovebox goes above the high pressure safety set point it will shut off.

NOTE: This is a dangerous situation and **caution** should be taken.




LC TECHNOLOGY SOLUTIONS INC.

2C Fanaras Drive, Salisbury, MA 01952 • Phone: 978-255-1620 • Fax: 978-428-0222
E-mail: info@lctechinc.com • Web: www.lctechinc.com

Section 22: Manual Purge

Purging with a Manual Purge Valve



LC TECHNOLOGY SOLUTIONS EL103


CHAMBER PRESSURE: +0000.0 mBAR
CHAMBER OXYGEN LEVEL: +0000.0 PPM
CHAMBER MOISTURE LEVEL: +0000.0 PPM

OFF
CHAMBER LIGHTS
OFF
CHAMBER PRESSURE CONTROL
CHAMBER CIRCULATION
CHAMBER SETTINGS

ALARMS
LOADLOCK1
PURGE
REGEN

F1 F2 F3 F4 F5 F6

1. From Main Screen verify Chamber Circulation is turned off. **OFF** should be displayed.



CHAMBER SETTINGS:

CHAMBER PRESSURE CONTROL HIGH SP: +0000.0 mBAR
CHAMBER PRESSURE CONTROL LOW SP: +0000.0 mBAR
CHAMBER GAS SOLENOID OFF DELAY SP: +00000 x.1SEC
CHAMBER VACUUM SOLENOID OFF DELAY SP: +00000 x.1SEC
CHAMBER PRESSURE HIGH ALARM SP: +0000.0 mBAR
CHAMBER PRESSURE LOW ALARM SP: +0000.0 mBAR
CHAMBER OXYGEN LEVEL HIGH ALARM SP: +0000.0 PPM
CHAMBER MOISTURE LEVEL HIGH ALARM SP: +0000.0 PPM

< RETURN CHAMBER SELECT CHAMBER 1 & 2

F1 F2 F3 F4 F5 F6

2. From Main Screen, press Chamber Settings to check pressure settings. Pressure settings should be set to positive values.

To verify settings for Purging follow these steps:

- a. Press Chamber Settings button.
- b. Make sure chamber pressure set points are set to **5 mBAR high and 2mBAR low.**



LC TECHNOLOGY SOLUTIONS INC.

2C Fanaras Drive, Salisbury, MA 01952 • Phone: 978-255-1620 • Fax: 978-428-0222
E-mail: info@lctechinc.com • Web: www.lctechinc.com



3. Open manual purge valve until you hear gas inlet valve come on and stay on.

4. For initial purging (glovebox is at room air) refer to Section 12 for the amount of time and gas.

5. When purging is complete follow these steps:
 - a. Close Purge Valve.
 - b. Return pressure set points to **2.5 mBAR high and 0.5 mBAR low.**
 - c. Start/Restart circulation by pressing Chamber Circulation button. **ON** should be displayed.



Section 23: Maintenance Schedule & Recommended Spare Parts

1. The gloves and glove O-rings should be changed once every three (3) months or as needed based on the condition of the gloves.
2. Vacuum pump oil should be changed at least every three (3) months and after a regeneration.
3. The inlet/outlet filters should be replaced every three (3) months.
4. The large antechamber door O-rings should be replaced as needed.
5. The small antechamber door O-rings should be replaced as needed.
6. If the system is equipped with a solvent removal system please follow the schedule below.
 - a. Internal solvent removal trap change every month; Part No. SR-101.
 - b. External solvent removal trap change charcoal every three (3) months; Part No. FM-018.

NOTE: With heavy solvent usage change the charcoal more frequently.

Spare Parts Listing

Part No.	Description	Price	Maintenance Schedule
GL-003	Gloves, Left/Right, 1 pair (Universal can be used in any Glovebox)	\$165.00/pair.	3 Months
OR-028	Glove O-Rings (4 per pair of gloves)	\$15.00 ea.	3 Months
SR-101	Internal Charcoal Trap (Set of 12)	\$450.00/set	Every Month
OR-111	Large Antechamber Door O-Ring (2 per door) (LC-1 Glovebox Systems)	\$25.00 ea.	As Needed
OR-110	Small Antechamber Door O-Ring (2 per door) (LC-1 Glovebox Systems)	\$10.00 ea.	As Needed
FL-102	Inlet/Outlet HEPA Filter	\$39.95 ea.	3 Months
FM-018	Activated Carbon for Solvent Removal System(10lbs per column)	\$10.00/lb	3 Months
FM-900	LC-1 Complete Filter Column Change	\$595.00/charge	As Needed
AN-023	Replacement Sensor for Oxygen Analyzer for Model OXY-IQ	\$300.00 ea.	As Needed




LC TECHNOLOGY SOLUTIONS INC.

2C Fanaras Drive, Salisbury, MA 01952 • Phone: 978-255-1620 • Fax: 978-428-0222
E-mail: info@lctechinc.com • Web: www.lctechinc.com

Section 24: Operational Accessories

24.1 Laminar Flow



1. From Main Screen press Laminar Flow Units button.




2. Laminar Flow Units Control screen will display.




LC TECHNOLOGY SOLUTIONS INC.

2C Fanaras Drive, Salisbury, MA 01952 • Phone: 978-255-1620 • Fax: 978-428-0222

E-mail: info@lctechinc.com • Web: www.lctechinc.com



3. Push Laminar Flow Unit button to turn laminar flow ON or OFF.




4. Press Return button to return to the main screen.



LC TECHNOLOGY SOLUTIONS INC.


2C Fanaras Drive, Salisbury, MA 01952 • Phone: 978-255-1620 • Fax: 978-428-0222
E-mail: info@lctechinc.com • Web: www.lctechinc.com

24.2 UV Ozone Cleaner



1. Open lid.

2. Place substrate in UV ozone cleaner.




3. Close lid and secure with knob.




LC TECHNOLOGY SOLUTIONS INC.

2C Fanaras Drive, Salisbury, MA 01952 • Phone: 978-255-1620 • Fax: 978-428-0222
E-mail: info@lctechinc.com • Web: www.lctechinc.com



4. From main screen press UV Ozone Cleaner button.


UV Ozone Cleaner Settings and Status screen will display.



5. Set desired cleaning and purging time by pressing the white box to right of set point.

The factory default for UV Cleaning Time SP: is +5.

The factory default for UV Chamber Purge Time SP: is +2.




6. Press green button on front of UV ozone cleaner to start process. Button will illuminate.



LC TECHNOLOGY SOLUTIONS INC.

2C Fanaras Drive, Salisbury, MA 01952 • Phone: 978-255-1620 • Fax: 978-428-0222
E-mail: info@lctechinc.com • Web: www.lctechinc.com

During the UV cleaning cycle the following information can be found on the UV Ozone Cleaner Settings And Status screen.



- A. UV Clean Cycle – This indicates if the cycle is ON or OFF.
- B. UV Clean Cycle Status – This is the current cycle; Cleaning or Purging.
- C. UV Cleaning Time Remaining – This is the amount of time remaining in the cleaning process. The time will tick down until it reaches zero.
- D. UV Chamber Purge Time Remaining – This is the amount of time remaining in the purge process. The time will tick down until it reaches zero.


Once the system reaches zeros (0) you have reached the end of the cycle and the process will automatically turn off.



LC TECHNOLOGY SOLUTIONS INC.

2C Fanaras Drive, Salisbury, MA 01952 • Phone: 978-255-1620 • Fax: 978-428-0222
E-mail: info@lctechinc.com • Web: www.lctechinc.com

24.3 Precision Hot Plate



1. Place substrate on the hot plate.



2. Turn power on by pushing switch down.



3. Press vacuum button to vacuum substrate down.



LC TECHNOLOGY SOLUTIONS INC.

2C Fanaras Drive, Salisbury, MA 01952 • Phone: 978-255-1620 • Fax: 978-428-0222
E-mail: info@lctechinc.com • Web: www.lctechinc.com



4. To adjust the temperature press and hold down the down arrow key.
5. SP1 will flash.
6. Set SP1 using the up and down arrow keys.
7. Once desired setting has been reached press P key.


NOTE: For advanced setting information refer to the vendor supplied manual located in the documentation package.




LC TECHNOLOGY SOLUTIONS INC.

2C Fanaras Drive, Salisbury, MA 01952 • Phone: 978-255-1620 • Fax: 978-428-0222
E-mail: info@lctechinc.com • Web: www.lctechinc.com

24.4 Spin Coater



1. Open lid.



2. Place substrate on spin coater chuck.

3. Close lid.

4. Press start to run.

NOTE: For detailed operating instructions refer to the vendor supplied manual located in the documentation package.





LC TECHNOLOGY SOLUTIONS INC.

2C Fanaras Drive, Salisbury, MA 01952 • Phone: 978-255-1620 • Fax: 978-428-0222
E-mail: info@lctechinc.com • Web: www.lctechinc.com

24.5 UV Press



1. Turn light source on.
2. Open UV press lid.




3. Load substrate into substrate holder.




LC TECHNOLOGY SOLUTIONS INC.


2C Fanaras Drive, Salisbury, MA 01952 • Phone: 978-255-1620 • Fax: 978-428-0222
E-mail: info@lctechinc.com • Web: www.lctechinc.com



4. Turn on suction.



5. Place cover glass in vacuum chuck.




6. Place vacuum chuck in press.



LC TECHNOLOGY SOLUTIONS INC.


2C Fanaras Drive, Salisbury, MA 01952 • Phone: 978-255-1620 • Fax: 978-428-0222
E-mail: info@lctechinc.com • Web: www.lctechinc.com



7. Close cover.



8. Press UV Press button on main screen.




9. UV Press Setting and Status screen will be displayed.



LC TECHNOLOGY SOLUTIONS INC.

2C Fanaras Drive, Salisbury, MA 01952 • Phone: 978-255-1620 • Fax: 978-428-0222
E-mail: info@lctechinc.com • Web: www.lctechinc.com



10. Set desired set point by pushing white box to the right of set point.


The factory default setting for the Chamber Evacuation Pressure SP is 750 mBAR.

The factory default setting for the UV Curing Time SP is 120 SEC.

The factory default setting for the Chamber Refill Time SP is 30 SEC.



11. Verify UV press clamping pressure is set at 30 psi.




12. Press green button on front of UV press to start process.



LC TECHNOLOGY SOLUTIONS INC.

2C Fanaras Drive, Salisbury, MA 01952 • Phone: 978-255-1620 • Fax: 978-428-0222
E-mail: info@lctechinc.com • Web: www.lctechinc.com

During the UV press cycle the following information can be found on the UV Press Settings And Status screen.



- A. UV Press Cycle – This indicates if the cycle is ON or OFF.
- B. UV Press Cycle Status – This is the current cycle; EVAC, Clamp, Curing, Refilling, Off.
- C. Chamber Pressure – This is the actual chamber pressure readout. Once the chamber reaches the set point it will change steps.
- D. UV Curing Time Remaining – This is the amount of time remaining in the curing process. The time will tick down until it reaches zero.
- E. Chamber Refill Time Remaining – This is the amount of time remaining in the refill process. The time will tick down until it reaches zero.

Once the system reaches zeros (0) you have reached the end of the cycle and the process will automatically turn off.



24.6 FUJIFILM Dimatix Materials Printer



1. To operate printer turn power on.


NOTE: For detailed operating instructions refer to the vendor supplied manual located in the documentation package.



LC TECHNOLOGY SOLUTIONS INC.

2C Fanaras Drive, Salisbury, MA 01952 • Phone: 978-255-1620 • Fax: 978-428-0222
E-mail: info@lctechinc.com • Web: www.lctechinc.com

24.7 Fisnar Glue Dispensing Robot



1. To operate the glue dispensing report turn power on.


NOTE: For detailed operating instructions refer to the vendor supplied manual located in the documentation package.



LC TECHNOLOGY SOLUTIONS INC.

2C Fanaras Drive, Salisbury, MA 01952 • Phone: 978-255-1620 • Fax: 978-428-0222
E-mail: info@lctechinc.com • Web: www.lctechinc.com


24.8 Heat & Cooling Plate





LC TECHNOLOGY SOLUTIONS INC.





2C Fanaras Drive, Salisbury, MA 01952 • Phone: 978-255-1620 • Fax: 978-428-0222
E-mail: info@lctechinc.com • Web: www.lctechinc.com



3. Turn on by pressing the switch to the ON position.
4. Turn off by pressing the switch to the OFF position.



To adjust the Set Point:

4. Press the  key on the temperature controller.
5. Use the  or  arrows to change set point.
6. Press the  key again to complete.

NOTE: Do not change any other settings on the temperature controller without first contacting LC Technology.

NOTE: The maximum operating temperature is 40°C to -100°C.



NOTE: Liquid nitrogen must be connected to the liquid nitrogen feed line for the plate to cool.

NOTE: Liquid nitrogen must not exceed 60 psi.



LC TECHNOLOGY SOLUTIONS INC.

2C Fanaras Drive, Salisbury, MA 01952 • Phone: 978-255-1620 • Fax: 978-428-0222
E-mail: info@lctechinc.com • Web: www.lctechinc.com

24.9 Heated Tray

1. Load items into the antechamber, for instructions see **Section 14**.



2. Once items are in the antechamber and the antechamber is evacuated, you can turn Heated Tray Control Power switch to the ON position.



3. Press the run process button on the oven controller.

See Oven Controller Setup Instructions below for more information.





LC TECHNOLOGY SOLUTIONS INC.

2C Fanaras Drive, Salisbury, MA 01952 • Phone: 978-255-1620 • Fax: 978-428-0222
E-mail: info@lctechinc.com • Web: www.lctechinc.com

4. Once process is complete refill the antechamber, prior to opening the door.

CAUTION! Contents will be hot.



LC TECHNOLOGY SOLUTIONS INC.

2C Fanaras Drive, Salisbury, MA 01952 • Phone: 978-255-1620 • Fax: 978-428-0222
E-mail: info@lctechinc.com • Web: www.lctechinc.com

24.9a Heated Tray Controller Set-Up

1. Turn power on.

NOTE: This is a ramp and hold programmable controller. It has nine (9) program segments. To program them follow the steps below.

2. Press scroll button.

Arrow Keys

Scroll

HEATED TRAY CONTROL



3. A green #1 will appear on the left side of this screen. Select desired temperature by using the ↑ and ↓ arrows. This is your first temperature set point.




4. Press scroll button.

5. Set desired time to reach required temperature. Use ↑ and ↓ arrows. This is your first ramp rate.

NOTE: When setting time hours are to the left of decimal point and minutes are to the right.

6. Press scroll button.

7. A green #2 will appear on this screen. You will once again select desired temperature using the ↑ and ↓ arrows. The same temperature should be used that was used in step 3. This is your first bake time.





LC TECHNOLOGY SOLUTIONS INC.

2C Fanaras Drive, Salisbury, MA 01952 • Phone: 978-255-1620 • Fax: 978-428-0222
E-mail: info@lctechinc.com • Web: www.lctechinc.com

8. Press scroll button. Select the desired amount of time you wish to bake items. Use ↑ and ↓ arrows.

9. Press Scroll Button.

10. There are nine (9) set values that can be utilized, if necessary. They follow the same pattern. If you were to use them the third set value would be a ramp up time and the fourth set value would be a bake or soak time.

11. If only set values 1 and 2 are needed, leave zeros (0) for the other set values. You will cycle through these set values by using the scroll button.

12. To save settings and return to main screen, press and hold the scroll button until a solid red number is displayed. This is the temperature inside the oven at the present moment.



13. You are now ready to run the oven.

NOTE: To cancel run process, press and hold run button until the green number on the left side of screen vanishes.





LC TECHNOLOGY SOLUTIONS INC.

2C Fanaras Drive, Salisbury, MA 01952 • Phone: 978-255-1620 • Fax: 978-428-0222
E-mail: info@lctechinc.com • Web: www.lctechinc.com


24.10 Solvent Sensor



If you purchased the optional solvent sensor your gas purifier will have a solvent sensor attached to it. (Pictured to the left)

The solvent sensor is designed to sense solvent vapor within the glovebox. The indication on the sensor is used to determine the condition of the activated carbon/ molecular sieve material used in your solvent removal system.

The sensor is sensitive to a variety of solvent vapor.



The content of solvent vapor is continuously displayed on the main plc screen.

The typical reading is 0.0.

If the reading is higher than 0.0 this is an indication that you need to change your activated carbon or reactivate the molecular sieve in the solvent removal filter column.


Note: We still recommend changing the filter material every 6 months regardless of the sensor reading.



LC TECHNOLOGY SOLUTIONS INC.

2C Fanaras Drive, Salisbury, MA 01952 • Phone: 978-255-1620 • Fax: 978-428-0222

E-mail: info@lctechinc.com • Web: www.lctechinc.com



By pressing the button on the main screen that says Solvent Monitor you will see the screen to the left. It displays the current solvent content inside the glovebox. It also displays the user adjustable set point for the solvent sensor alarm. This is typically set to 0.1 as you will want to get an alarm as soon as there is excess solvent vapor inside the glovebox.


If the solvent content is above the alarm set point it will trigger an alarm on the alarm screen. This is an indication that you need to replace or reactivate the solvent removal filter column material.



LC TECHNOLOGY SOLUTIONS INC.

2C Fanaras Drive, Salisbury, MA 01952 • Phone: 978-255-1620 • Fax: 978-428-0222
E-mail: info@lctechinc.com • Web: www.lctechinc.com


24.11 Power Saver Mode



The power saver mode is designed to switch off the vacuum pump when it is not needed to run the system. This saves both energy and cuts down on noise.

To turn on power saver mode press the display button labeled Power Save Mode.

This will take you to the below screen.



Once in this screen you should press the Power Save Mode to start this feature. ON will be displayed

The pump should turn off, (this is normal). Then you should enter the time that you want the pump to run for (off delay time). Typically this would be set to 60 Minutes.

The system will continue to monitor the pressure in the glovebox to determine if the pump needs to be activated.

The pump should turn off automatically once it has reached the preset vacuum level.



LC TECHNOLOGY SOLUTIONS INC.

2C Fanaras Drive, Salisbury, MA 01952 • Phone: 978-255-1620 • Fax: 978-428-0222
E-mail: info@lctechinc.com • Web: www.lctechinc.com

Section 25: Technical Support Contact Information from LC Technology Solutions & Partners

<p>LC TECHNOLOGY SOLUTIONS INC.</p> <p>Tel: (978) 255-1620 Cell: (978) 992-1729 or 978-289-7723 Fax: (978) 428-0222 Email: info@lctechinc.com Web: www.lctechinc.com</p> <p><u>It helps to provide your Project #, which can be found on the back of your glovebox.</u></p> <p>AGILENT (Scroll & Vacuum Pumps)</p> <p>Alan Bird Tel: (978) 387-3035 Email: alan.bird@agilent.com</p> <p>Ray Hula Tel: (800) 882-7426, #3 for Tech Support Email: ray.hula@agilent.com</p> <p>James Ramsden Tel: (800) 882-7426, #3 for Tech Support Email: james.ramsden@agilent.com Web: www.agilent.com</p> <p>EDWARDS (RV3, 12 Vacuum & Scroll Pumps)</p> <p>Randy Morse Tel: (800) 848-9800 Ext. 3459</p> <p>General Tech Support Tel: (800) 848-9800 Ext. 3344</p> <p>FISNAR (Robotic Components)</p> <p>Shailesh Lad Phone: (973) 646-5044 Ext. 1302 Email: slad@fisnar.com Web: www.fisnar.com</p>	<p>FUJIFILM DIMATIX, INC. (Printers & Robotic Components)</p> <p>Stephanie Scattareggia or Harrison Nguyen Tel: (408) 565-7025 Fax: (408) 565-7060 Email: sscattareggia@fujifilm.com Email: hnguyen@fujifilm.com Web: www.fujifilm.com</p> <p>GARDNER DENVER (Welch Vacuum) (Solvent Purifiers)</p> <p>Frank Dziedzic Tel: (847) 588-2365 Email: frank.dziedzic@gardnerdenver.com</p> <p>Mark Suda Tel: (847) 588-2358 Email: Mark.Suda@gardnerdenver.com</p> <p>GE SENSING (Moisture & Oxygen Analyzers)</p> <p>Tel: (800) 833-9438 Email: sensing@ge.com Web: www.ge.com</p> <p>SCS - Specialty Coating Systems (Spin Coaters & Hot Plates)</p> <p>Shawn Gordon Tel: (317) 472-1223 Email: sgordon@scscoatings.com Web: www.scscoatings.com</p> <p>TROVATO (Evaporation Systems)</p> <p>Tom Trovato Tel: (585) 742-8070 Fax: (585) 742-3811 Email: tom@trovato.org Web: www.trovato.org</p>
--	---