



LC TECHNOLOGY SOLUTIONS INC.

2C Fanaras Drive, Salisbury, MA 01952 • Phone: 978-255-1620 • Fax: 978-428-0222
E-mail: info@lctechinc.com • Web: www.lctechinc.com

SP-1 Stand Alone Solvent Purification System

Operation Manual





LC TECHNOLOGY SOLUTIONS INC.

2C Fanaras Drive, Salisbury, MA 01952 • Phone: 978-255-1620 • Fax: 978-428-0222
E-mail: info@lctechinc.com • Web: www.lctechinc.com


Table of Contents

Section 1: System Overview	3
Section 2: Ordering Solvents.....	5
Section 3: Installation Instructions.....	6
3.1 System Piping Connections.....	6
3.2 Special Instructions for Methanol.....	7
3.3 Special Instructions for THF, Ether & Dioxane.....	9
3.4 Utility Connections.....	10
3.5 Filling Solvent Kegs.....	12
Section 4: Bubble Degas Procedure	14
Section 5: Operational Instructions.....	16
5.1 Initial Start-Up.....	16
5.2 Operation of Solvent Purification System.....	17
Section 6: Technical Support Contact Information from LC Technology Solutions & Partners...	20



SP-1 Stand Alone Solvent Purification Operation Manual

Section 1: System Overview






LC TECHNOLOGY SOLUTIONS INC.

2C Fanaras Drive, Salisbury, MA 01952 • Phone: 978-255-1620 • Fax: 978-428-0222
E-mail: info@lctechinc.com • Web: www.lctechinc.com

Solvent
Connection



Gas Connection

Pressure Relief Valve

The Project Number is located on a sticker on the back of your system. Please have this number available when calling for service information.



Section 2: Ordering Solvents

Solvents need to be ACS grade with a starting water content of 250 ppm or less. If you are using THF or ether it must be inhibitor free.

The following list includes some vendors who provide these types of solvents.

Sigma Aldrich

http://www.sigmaaldrich.com/Area_of_Interest/Research_Essentials/Solvents.html

Avantor Materials

<http://www.avantormaterials.com/>


Fisher Scientific

<http://www.fishersci.com/>



Section 3: Installation Instructions

3.1 System Piping Connections



1. Connect all inter-connecting piping supplied with system. Connections will be labeled by the factory.
2. Match the corresponding numbers on component with the corresponding number on system.


NOTE: See system pictures for connection details.




3.2 Special Instructions for Methanol

NOTE! Your system has been marked with a Methanol designated channel and solvent keg. You must use this channel only for Methanol.


If you have not told us that you will be using Methanol in your system, please contact LC Technology. We will provide you with the appropriate filter for use with your system.



1. After filling solvent keg with Methanol remove right side panel.



2. Locate the Methanol Filter.



3. Gently twist open small white bleed valve on top of filter.
NOTE: Do not remove cap completely!
Methanol level will rise and fill entire filter.
Once filled retighten bleed valve.

4. Reattach right side panel.

5. Resume normal operation.



LC TECHNOLOGY SOLUTIONS INC.

2C Fanaras Drive, Salisbury, MA 01952 • Phone: 978-255-1620 • Fax: 978-428-0222
E-mail: info@lctechinc.com • Web: www.lctechinc.com

NOTE: It is very important to fill the system slowly with methanol. Methanol has a tendency to cause the filter column to heat up. This excessive heat can cause the filter material to breakdown.

To fill the system slowly, open the dispensing valve for about 30 seconds and then close the valve. Check the filter columns to see if they are warm. If they are not warm, open the dispensing valve for another 30 seconds and check again. If they are warm wait for them to cool and then dispense again. Do this until the methanol begins to flow out of the dispensing port.



LC TECHNOLOGY SOLUTIONS INC.

2C Fanaras Drive, Salisbury, MA 01952 • Phone: 978-255-1620 • Fax: 978-428-0222
E-mail: info@lctechinc.com • Web: www.lctechinc.com

3.3 Special Instructions for THF, Ether & Dioxane

NOTE: You must use inhibitor-free solvent when using THF, Ether and/or Dioxane.



LC TECHNOLOGY SOLUTIONS INC.

2C Fanaras Drive, Salisbury, MA 01952 • Phone: 978-255-1620 • Fax: 978-428-0222
E-mail: info@lctechinc.com • Web: www.lctechinc.com

3.4 Utility Connections



Inert Gas Connection

1. Connect inert gas supply to the Gas Inlet connection.
 - 3/8" Hose barb included.
 - Recommend connecting to 3/8" reinforced tygon tubing.
 - Inert gas supply should be nitrogen or argon at 15 psi.

Power Connection for the Vacuum Pump

US locations plug the vacuum pump into 115V power supply.

International locations plug the vacuum pump into 220V power supply.



Important Safety Note

Verify the system is grounded properly.




Connect Vacuum Pump Exhaust to Fume Hood

- 3/8" Hose barb included.
- Recommend connecting to 3/8" reinforced tygon tubing.




LC TECHNOLOGY SOLUTIONS INC.

2C Fanaras Drive, Salisbury, MA 01952 • Phone: 978-255-1620 • Fax: 978-428-0222
E-mail: info@lctechinc.com • Web: www.lctechinc.com



Connect Fume Hood Connection

1. Connect the fume hood connection located on top of the solvent purifier to the fume hood.
 - This connection is optional.
 - 4" Hole provided.



Locking Nut

Set Screw

Setting Gas Regulator / Pressure

The inlet gas regulator has been set by the factory at 7 psi.

If you need to adjust this, adjust by turning the set screw on the regulator.

NOTE: Loosen the locking nut before adjusting set screw.



3.5 Filling Solvent Kegs



1. Remove solvent keg from fire cabinet.



2. Disconnect gas and solvent lines.
Lines are equipped with quick connects.
3. Remove by pulling gently back on the ridged portion of the quick connect.

3. Take solvent keg to fume hood or other well-ventilated area for filling.



4. Loosen clamp at top of keg.
NOTE: Relieve gas pressure in keg prior to removing.
5. Remove clamp.



6. Lift and remove solvent connection flange and O-ring.



LC TECHNOLOGY SOLUTIONS INC.

2C Fanaras Drive, Salisbury, MA 01952 • Phone: 978-255-1620 • Fax: 978-428-0222
E-mail: info@lctechinc.com • Web: www.lctechinc.com



7. Pour solvent into keg.



8. Replace the O-ring and solvent connection flange onto the solvent keg.
9. Replace and tighten the clamp.

10. Reattach gas connection and then the solvent connection.

NOTE: Be sure the color coded labels attached to the gas and solvent lines match the color coded label affixed to the front of the solvent keg.



11. Place the solvent keg back in the fire cabinet.

12. Repeat the above steps as necessary for filling additional solvent kegs.

NOTE: Once you have filled the solvent purifier with a particular type of solvent you cannot change to a different solvent without changing the filter columns and completely cleaning the entire system. If you do not change the filters and clean the system you will experience cross contamination of the solvents. Please contact LC Technology if you decide to change from one solvent to another.




Section 4: Bubble Degas Procedure

1. Bubble degassing should be performed in a fume hood.


2. Disconnect solvent keg from system and place in it in fume hood.

3. Connect 1/4" Swagelok fitting of female quick connect to an inert gas supply.



Male Quick Connect

4. Connect female quick connect to male quick connection on the solvent keg. This will supply inert gas down the dip tube to the bottom of the solvent keg. (Inert Gas Line)





LC TECHNOLOGY SOLUTIONS INC.

2C Fanaras Drive, Salisbury, MA 01952 • Phone: 978-255-1620 • Fax: 978-428-0222
E-mail: info@lctechinc.com • Web: www.lctechinc.com



Female Quick Connect

5. Connect male quick connect to the female quick connect on the solvent keg. 1/4" Swagelok fitting is provided to connect a vent line, if necessary. Solvent vapor will be coming out this connection and should be vented. (Vent Line)



6. Turn on inert gas supply to approximately (5) five psi. This will force inert gas down the dip tube and let it bubble up through the solvent and out the vent line.

7. Bubble degas for (5) five minutes or desired amount of time.

8. Turn off inert gas supply.

9. Disconnect bubble degas lines and remove from solvent keg.

10. Reattach solvent keg to system. Make sure all lines match the appropriate color coded connections.



Section 5: Operational Instructions

5.1 Initial Start-Up

1. After assembly of system, filling solvent keg(s) and bubble degassing turn on vacuum pump and inert gas supply.



2. Take an empty beaker, flask, or bottle and place it beneath the 24/40 connector.



3. Turn the lower operational valve to Dispense.

4. Leave lower operational valve in dispense position until a continuous flow of solvent is attained.

NOTE: First gas will empty quickly, then slow and finally the solvent will flow. This process can take up to 1-2 minutes.




5.2 Operation of Solvent Purification System




1. Turn on the vacuum pump and inert gas supply.


2. Select which solvent to collect.



3. Make sure the operational valves are in the Closed position.




4. Place jack stand, cork ring and collection vessel beneath the appropriate 24/40 connector.




5. Raise the jack stand until the collection vessel connects to the 24/40 connector.

NOTE: Make sure the connection is tight.



6. Open valve on top of the collection vessel.



7. Close valve on the syringe port.



8. Turn the lower operational valve to Vac/Gas.



LC TECHNOLOGY SOLUTIONS INC.

2C Fanaras Drive, Salisbury, MA 01952 • Phone: 978-255-1620 • Fax: 978-428-0222
E-mail: info@lctechinc.com • Web: www.lctechinc.com



9. Turn the upper operational valve to Vacuum.

NOTE: The vacuum level will be displayed on the vacuum gauge. The recommended vacuum level should be -28.5 in. Hg.



10. Once correct level of vacuum is reached turn the upper operational valve to Gas for one second then return to Closed.

11. **Repeat steps 9-10 two more times for a total of three cycles.



12. Turn upper operational valve to vacuum.

NOTE: Evacuate collection vessel to -28.5 in. Hg.



13. Turn upper operational valve to Closed.



14. Turn lower operational valve to dispense.

15. Collect the desired amount of solvent.



16. Turn lower operational valve to Vac/Gas position.



LC TECHNOLOGY SOLUTIONS INC.

2C Fanaras Drive, Salisbury, MA 01952 • Phone: 978-255-1620 • Fax: 978-428-0222
E-mail: info@lctechinc.com • Web: www.lctechinc.com



17. Turn upper operational valve to Gas for several seconds.



18. Turn upper and lower operational valves to Closed.



19. Close upper valve on collection vessel.

20. Remove collection vessel by lowering the jack stand.



Section 6: Technical Support Contact Information from LC Technology Solutions & Partners

<p>LC TECHNOLOGY SOLUTIONS INC.</p> <p>Tyler Kapsimalis Tel: (978) 255-1620 Ext. 107 Cell: (978) 992-1729 Fax: (978) 428-0222 Email: tkapsimalis@lctechinc.com Web: www.lctechinc.com</p> <p>It helps to provide your Project #, which can be found on the back of your glovebox.</p> <p>AGILENT (Scroll & Vacuum Pumps)</p> <p>Ray Hula Tel: (800) 882-7426, #3 for Tech Support Email: ray.hula@agilent.com</p> <p>James Ramsden Tel: (800) 882-7426, #3 for Tech Support Email: james.ramsden@agilent.com Web: www.agilent.com</p> <p>EDWARDS (RV3 & 12 Vacuum Pumps)</p> <p>Randy Morse Tel: (800) 848-9800 Ext. 3459</p> <p>General Tech Support Tel: (800) 848-9800 Ext. 3344</p> <p>FISNAR (Robotic Components)</p> <p>Shailesh Lad Phone: (973) 646-5044 Ext. 1302 Email: slad@fisnar.com Web: www.fisnar.com</p>	<p>FUJIFILM DIMATIX, INC. (Printers & Robotic Components)</p> <p>Stephanie Scattareggia or Harrison Nguyen Tel: (408) 565-7025 Fax: (408) 565-7060 Email: sscattareggia@fujifilm.com Email: hnguyen@fujifilm.com Web: www.fujifilm.com</p> <p>GARDNER DENVER (Welch Vacuum) (Solvent Purifiers)</p> <p>Frank Dziedzic Tel: (847) 588-2365 Email: frank.dziedzic@gardnerdenver.com</p> <p>Mark Suda Tel: (847) 588-2358 Email: Mark.Suda@gardnerdenver.com</p> <p>GE SENSING (Moisture & Oxygen Analyzers)</p> <p>Tel: (800) 833-9438 Email: sensing@ge.com Web: www.ge.com</p> <p>SCS - Specialty Coating Systems (Spin Coaters & Hot Plates)</p> <p>Shawn Gordon Tel: (317) 472-1223 Email: sgordon@scscoatings.com Web: www.scscoatings.com</p> <p>TROVATO (Evaporation Systems)</p> <p>Tom Trovato Tel: (585) 742-8070 Fax: (585) 742-3811 Email: tom@trovato.org Web: www.trovato.org</p>
--	---