



LC TECHNOLOGY SOLUTIONS INC.

2C Fanaras Drive, Salisbury, MA 01952 • Phone: 978-255-1620 • Fax: 978-428-0222  
E-mail: info@lctechinc.com • Web: www.lctechinc.com

## SPBT-1 Bench Top Solvent Purification System

### Operation Manual





# LC TECHNOLOGY SOLUTIONS INC.

2C Fanaras Drive, Salisbury, MA 01952 • Phone: 978-255-1620 • Fax: 978-428-0222  
E-mail: [info@lctechinc.com](mailto:info@lctechinc.com) • Web: [www.lctechinc.com](http://www.lctechinc.com)


## Table of Contents

Section 1: System Overview .....	3
Section 2: Ordering Solvents.....	6
Section 3: Installation Instructions .....	7
3.1 System Piping Connections .....	7
3.2 Special Instructions for Methanol.....	8
3.3 Special Instructions for THF & Ether.....	9
3.4 Utility Connections.....	10
3.5 Filling Solvent Kegs .....	12
Section 4: Bubble Degas Procedure .....	15
Section 5: Operational Instructions.....	17
5.1 Initial Start-Up.....	17
5.2 Operation of Solvent Purification System.....	18
Section 6: Technical Support Contact Information from LC Technology Solutions & Partners...	21



# SPBT-1 Bench Top Solvent Purification Operation Manual

## Section 1: System Overview





# LC TECHNOLOGY SOLUTIONS INC.

2C Fanaras Drive, Salisbury, MA 01952 • Phone: 978-255-1620 • Fax: 978-428-0222  
E-mail: [info@lctechinc.com](mailto:info@lctechinc.com) • Web: [www.lctechinc.com](http://www.lctechinc.com)




Gas Regulator



# LC TECHNOLOGY SOLUTIONS INC.

2C Fanaras Drive, Salisbury, MA 01952 • Phone: 978-255-1620 • Fax: 978-428-0222  
E-mail: [info@lctechinc.com](mailto:info@lctechinc.com) • Web: [www.lctechinc.com](http://www.lctechinc.com)



**The Project Number is located on a sticker on the back of your system. Please have this number available when calling for service information.**



---

## Section 2: Ordering Solvents

---

Solvents need to be ACS grade with a starting water content of 250 ppm or less. If you are using THF or ether it must be inhibitor free.

The following list includes some vendors who provide these types of solvents.

Sigma Aldrich

[http://www.sigmaaldrich.com/Area\\_of\\_Interest/Research\\_Essentials/Solvents.html](http://www.sigmaaldrich.com/Area_of_Interest/Research_Essentials/Solvents.html)

Avantor Materials

<http://www.avantormaterials.com/>


Fisher Scientific

<http://www.fishersci.com/>



## Section 3: Installation Instructions

### 3.1 System Piping Connections



1. Connect all inter-connecting piping supplied with system. Connections will be labeled by the factory.
2. Match the corresponding numbers on component with the corresponding number on system.


**NOTE:** See system pictures for connection details.




### 3.2 Special Instructions for Methanol

**NOTE!** Your system has been marked with a Methanol designated channel and solvent keg. You must use this channel only for Methanol.


If you have not told us that you will be using Methanol in your system, please contact LC Technology. We will provide you with the appropriate filter for use with your system.



1. After filling solvent keg with Methanol remove right side panel.



2. Locate the Methanol Filter.



3. Gently twist open small white bleed valve on top of filter.  
**NOTE: Do not remove cap completely!**  
Methanol level will rise and fill entire filter.  
Once filled retighten bleed valve.

4. Reattached right side panel.

5. Resume normal operation.





## LC TECHNOLOGY SOLUTIONS INC.


2C Fanaras Drive, Salisbury, MA 01952 • Phone: 978-255-1620 • Fax: 978-428-0222  
E-mail: [info@lctechinc.com](mailto:info@lctechinc.com) • Web: [www.lctechinc.com](http://www.lctechinc.com)

### 3.3 Special Instructions for THF & Ether

**NOTE:** You must use inhibitor-free solvent when using THF and/or Ether.



### 3.4 Utility Connections




#### Inert Gas Connection

1. Connect inert gas supply to the Gas Inlet connection.
  - 3/8" Hose barb included.
  - Recommend connecting to 3/8" reinforced tygon tubing.
  - Inert gas supply should be nitrogen or argon at 15 psi.

#### Power Connection for the Vacuum Pump


US locations plug the vacuum pump into 115V power supply.

International locations plug the vacuum pump into 220V power supply.



#### Important Safety Note

Verify the system is grounded properly.




#### Connect Vacuum Pump Exhaust to Fume Hood

- 3/8" Hose barb included.
- Recommend connecting to 3/8" reinforced tygon tubing.



## LC TECHNOLOGY SOLUTIONS INC.

2C Fanaras Drive, Salisbury, MA 01952 • Phone: 978-255-1620 • Fax: 978-428-0222  
E-mail: [info@lctechinc.com](mailto:info@lctechinc.com) • Web: [www.lctechinc.com](http://www.lctechinc.com)



Locking Nut

Set Screw

### Setting Gas Regulator / Pressure

The inlet gas regulator has been set by the factory at 7 psi.

If you need to adjust this, adjust by turning the set screw on the regulator.

**NOTE:** Loosen the locking nut before adjusting set screw.



### 3.5 Filling Solvent Kegs



1. Remove solvent keg from the top of the system.



2. Disconnect gas and solvent lines.
3. Lines are equipped with quick connects. Remove by pulling gently back on the ridged portion of the quick connect.

4. Take solvent keg to fume hood or other well-ventilated area for filling.



5. Loosen clamp at top of keg.  
**NOTE:** Relieve gas pressure in keg prior to removing.
6. Remove clamp.



7. Lift and remove solvent connection flange and O-ring.



## LC TECHNOLOGY SOLUTIONS INC.

2C Fanaras Drive, Salisbury, MA 01952 • Phone: 978-255-1620 • Fax: 978-428-0222  
E-mail: [info@lctechinc.com](mailto:info@lctechinc.com) • Web: [www.lctechinc.com](http://www.lctechinc.com)



1. Pour solvent into keg.



2. Replace the O-ring and solvent connection flange onto the solvent keg.
3. Replace and tighten the clamp.

4. Reattach gas connection and then the solvent connection.

**NOTE:** Be sure the color coded labels attached to the gas and solvent lines match the color coded label affixed to the front of the solvent keg.



5. Place the solvent keg back on top of the system.

6. Repeat the above steps as necessary for filling additional solvent kegs.

**NOTE:** Once you have filled the solvent purifier with a particular type of solvent you cannot change to a different solvent without changing the filter columns and completely cleaning the entire system. If you do not change the filters and clean the system you will experience cross



## LC TECHNOLOGY SOLUTIONS INC.

2C Fanaras Drive, Salisbury, MA 01952 • Phone: 978-255-1620 • Fax: 978-428-0222  
E-mail: [info@lctechinc.com](mailto:info@lctechinc.com) • Web: [www.lctechinc.com](http://www.lctechinc.com)

contamination of the solvents. Please contact LC Technology if you decide to change from one solvent to another.




## Section 4: Bubble Degas Procedure

1. Bubble degassing should be performed in a fume hood.


2. Disconnect solvent keg from system and place in it in fume hood.

3. Connect ¼" Swagelok fitting of female quick connect to an inert gas supply.




Male Quick Connect

4. Connect female quick connect to male quick connection on the solvent keg. This will supply inert gas down the dip tube to the bottom of the solvent keg. (Inert Gas Line)



Female Quick Connect


5. Connect male quick connect to the female quick connect on the solvent keg. ¼" Swagelok fitting is provided to connect a vent line, if nessassary. Solvent vapor will be coming out this connection and should be vented. (Vent Line)





## LC TECHNOLOGY SOLUTIONS INC.

2C Fanaras Drive, Salisbury, MA 01952 • Phone: 978-255-1620 • Fax: 978-428-0222  
E-mail: [info@lctechinc.com](mailto:info@lctechinc.com) • Web: [www.lctechinc.com](http://www.lctechinc.com)



6. Turn on inert gas supply to approximately (5) five psi. This will force inert gas down the dip tube and let it bubble up through the solvent and out the vent line.

7. Bubble degas for (5) five minutes or desired amount of time.

8. Turn off inert gas supply.

9. Disconnect bubble degas lines and remove from solvent keg.

10. Reattach solvent keg to system. Make sure all lines match the appropriate color coded connections.






## Section 5: Operational Instructions

### 5.1 Initial Start-Up

1. After assembly of system, filling solvent keg(s) and bubble degassing turn on vacuum pump and inert gas supply.



2. Take an empty beaker, flask, or bottle and place it beneath the 24/40 connector.



3. Turn the lower operational valve to Dispense.

4. Leave lower operational valve in dispense position until a continuous flow of solvent is attained.

**NOTE:** First gas will empty quickly, then slow and finally the solvent will flow. This process can take up to 1-2 minutes.








5. Solvent keg supplied with system hold five (5) liters of solvent. After initial start-up you may have to refill the keg.

This should only occur at initial start-up.

**NOTE:** if you use Methanol in your system you may experience some cloudiness in your purified solvent. This is a small amount of dust on the filter material. You will need to flush (2) or (3) liters of Methanol through the system. This should clear up any remaining cloudiness in your Methanol. If it continues to be cloudy, please continue to flush additional Methanol through the system until it becomes clear.



## 5.2 Operation of Solvent Purification System

		<p>1. Turn on the vacuum pump and inert gas supply.</p>
		<p>2. Select which solvent to collect.</p>
		<p>3. Make sure the operational valves are in the Closed position.</p>
		<p>4. Place jack stand, cork ring and collection vessel beneath the appropriate 24/40 connector.</p>
		<p>5. Raise the jack stand until the collection vessel connects to the 24/40 connector.</p> <p><b>NOTE:</b> Make sure the connection is tight.</p>
		<p>6. Open valve on top of the collection vessel.</p>
		<p>7. Close valve on the syringe port.</p>
		<p>8. Turn the lower operational valve to Vac/Gas.</p>



## LC TECHNOLOGY SOLUTIONS INC.

2C Fanaras Drive, Salisbury, MA 01952 • Phone: 978-255-1620 • Fax: 978-428-0222  
E-mail: info@lctechinc.com • Web: www.lctechinc.com



9. Turn the upper operational valve to Vacuum.

**NOTE:** The vacuum level will be displayed on the vacuum gauge. The recommended vacuum level should be - 28.5 in. Hg.



10. Once correct level of vacuum is reached turn the upper operational valve to Gas for one second.

11. \*\*Repeat steps 9-10 two more times for a total of three cycles.



12. Turn upper operational valve to vacuum.

**NOTE:** Evacuate collection vessel to - 28.5 in. Hg.



13. Turn upper operational valve to Closed.



14. Turn lower operational valve to dispense.

15. Collect the desired amount of solvent.



16. Turn lower operational valve to Vac/Gas position.



## LC TECHNOLOGY SOLUTIONS INC.

2C Fanaras Drive, Salisbury, MA 01952 • Phone: 978-255-1620 • Fax: 978-428-0222  
E-mail: [info@lctechinc.com](mailto:info@lctechinc.com) • Web: [www.lctechinc.com](http://www.lctechinc.com)



17. Turn upper operational valve to Gas for several seconds.



18. Turn upper and lower operational valves to Closed.



19. Close upper valve on collection vessel.

20. Remove collection vessel by lowering the jack stand.



**Section 6: Technical Support Contact Information from LC Technology Solutions & Partners**

<p><b>LC TECHNOLOGY SOLUTIONS INC.</b></p> <p>Tyler Kapsimalis <b>Tel:</b> (978) 255-1620 Ext. 107 <b>Cell:</b> (978) 992-1729 <b>Fax:</b> (978) 428-0222 <b>Email:</b> <a href="mailto:tkapsimalis@lctechinc.com">tkapsimalis@lctechinc.com</a> <b>Web:</b> <a href="http://www.lctechinc.com">www.lctechinc.com</a></p> <p><a href="#">It helps to provide your Project #, which can be found on the back of your glovebox.</a></p> <p><b>AGILENT</b> <b>(Scroll &amp; Vacuum Pumps)</b></p> <p>Ray Hula <b>Tel:</b> (800) 882-7426, #3 for Tech Support <b>Email:</b> <a href="mailto:ray.hula@agilent.com">ray.hula@agilent.com</a></p> <p>James Ramsden <b>Tel:</b> (800) 882-7426, #3 for Tech Support <b>Email:</b> <a href="mailto:james.ramsden@agilent.com">james.ramsden@agilent.com</a> <b>Web:</b> <a href="http://www.agilent.com">www.agilent.com</a></p> <p><b>EDWARDS</b> <b>(RV3 &amp; 12 Vacuum Pumps)</b></p> <p>Randy Morse <b>Tel:</b> (800) 848-9800 Ext. 3459</p> <p>General Tech Support <b>Tel:</b> (800) 848-9800 Ext. 3344</p> <p><b>FISNAR</b> <b>(Robotic Components)</b></p> <p>Shailesh Lad <b>Phone:</b> (973) 646-5044 Ext. 1302 <b>Email:</b> <a href="mailto:slad@fisnar.com">slad@fisnar.com</a> <b>Web:</b> <a href="http://www.fisnar.com">www.fisnar.com</a></p>	<p><b>FUJIFILM DIMATIX, INC.</b> <b>(Printers &amp; Robotic Components)</b></p> <p>Stephanie Scattareggia or Harrison Nguyen <b>Tel:</b> (408) 565-7025 <b>Fax:</b> (408) 565-7060 <b>Email:</b> <a href="mailto:sscattareggia@fujifilm.com">sscattareggia@fujifilm.com</a> <b>Email:</b> <a href="mailto:hnguyen@fujifilm.com">hnguyen@fujifilm.com</a> <b>Web:</b> <a href="http://www.fujifilm.com">www.fujifilm.com</a></p> <p><b>GARDNER DENVER (Welch Vacuum)</b> <b>(Solvent Purifiers)</b></p> <p>Frank Dziedzic <b>Tel:</b> (847) 588-2365 <b>Email:</b> <a href="mailto:frank.dziedzic@gardnerdenver.com">frank.dziedzic@gardnerdenver.com</a></p> <p>Mark Suda <b>Tel:</b> (847) 588-2358 <b>Email:</b> <a href="mailto:Mark.Suda@gardnerdenver.com">Mark.Suda@gardnerdenver.com</a></p> <p><b>GE SENSING</b> <b>(Moisture &amp; Oxygen Analyzers)</b></p> <p><b>Tel:</b> (800) 833-9438 <b>Email:</b> <a href="mailto:sensing@ge.com">sensing@ge.com</a> <b>Web:</b> <a href="http://www.ge.com">www.ge.com</a></p> <p><b>SCS - Specialty Coating Systems</b> <b>(Spin Coaters &amp; Hot Plates)</b></p> <p>Shawn Gordon <b>Tel:</b> (317) 472-1223 <b>Email:</b> <a href="mailto:sgordon@scscoatings.com">sgordon@scscoatings.com</a> <b>Web:</b> <a href="http://www.scscoatings.com">www.scscoatings.com</a></p> <p><b>TROVATO</b> <b>(Evaporation Systems)</b></p> <p>Tom Trovato <b>Tel:</b> (585) 742-8070 <b>Fax:</b> (585) 742-3811 <b>Email:</b> <a href="mailto:tom@trovato.org">tom@trovato.org</a> <b>Web:</b> <a href="http://www.trovato.org">www.trovato.org</a></p>
--	---